For Research Use Only. Not for use in diagnostic procedures.

# **His tagged Protein PURIFICATION KIT**

(MoAb. clone 2D8)

CODE No. 3310

*PURIFICATION to Maintains Protein Activity  
from eukaryote cell lysate and culture supernatant*

## **Product Description**

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more. Thus, choosing a method of purification is an important aspect in maintaining protein structure and function.

The hexa Histidine-tag (6xHis tag) is one of the most common tags used to facilitate the purification of recombinant proteins. Metal chelate affinity chromatography is widely used for purification of His-tagged proteins. Unfortunately, some serum components are absorbed by the metal chelate affinity columns, making these columns impractical for the purification of proteins secreted into the culture supernatant. Indeed many proteins have intrinsic histidine residues or other untagged host proteins bind the metal chelate affinity columns and elute with His tagged protein. These untagged contaminants may be removed using an additional purification or laborious optimizing the imidazole concentration.

MBL's His tagged Protein PURIFICATION KIT is designed for the isolation of His tagged protein from cell culture supernatants containing serum and cell lysate under neutral pH condition. Severe conditions such as acidic or alkaline elution denature protein structure. However, a neutral pH elution can preserve protein activity and native conformation. MBL has developed the Anti-His tag Beads to purify His tagged proteins quickly and efficiently. As the Beads can be used at neutral pH, the purified proteins can maintain the activity and conformation. The elution of His tagged proteins from the Beads is achieved by the addition of the 6xHis peptide. As the 6xHis peptide competes with His tagged proteins on the Beads, the purified proteins do not lose the protein activity. The simple procedures of this kit have been optimized by using a Spin Column resulting in high efficiency.

## **Kit Components**

Components sufficient for conducting 20 times purifications of His tagged protein.

- |                              |   |
|------------------------------|---|
| <b>1. Anti-His tag Beads</b> | 25% slurry: 100 $\mu$ L beads in 400 $\mu$ L total volume in PBS with 0.09 % sodium azide as preservative |
| <b>2. Elution Peptide</b>    | 6xHis peptide, 2 mg in 1 mL PBS after reconstitution  |
| <b>3. Spin Columns Sets</b>  | 20 columns with pre-inserted bottom plugs and top caps  |
| <b>4. Wash Concentrate</b>   | 10x concentrate, 6 mL x 2 tubes   |

## **Storage**

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

## **Product Capacity**

The purification capacity of the Anti-His tag Beads varies depending upon the His tagged protein. For examples, 5  $\mu$ L of Anti-His Beads (20  $\mu$ L slurry) bound 5  $\mu$ g of a His tagged protein (32 kDa) and eluted 3  $\mu$ g of purified protein.

## **Materials Required but not Provided**

1. Microcentrifuge capable of 15,000 x g
2. Sampling tube (1.5 mL)
3. End-over-end rotator
4. PBS
5. Lysis buffer

Suitable Lysis buffer varies with cell type.

**Note:** see Additional Information

Homemade Lysis buffer

20-50 mM	Tris-HCl, pH 7.5 or HEPES-KOH, pH 7.5
50-250 mM	NaCl
5 mM	EDTA
1%	NP-40 or Triton X-100

if necessary add Protease Inhibitor Cocktail  
(e.g. SIGMA code P8340, PIERCE code 78415).

## **Protocols**

### **Introduction**

The kit is optimized under the native conditions only, and it is not recommended under denaturing conditions and also for purification of aggregated, unstable, and insoluble protein (e.g. inclusion bodies). Proteins solubilized with such as 6 M Guanidine-HCl or 8 M Urea can not be purified using this Kit (see Additional Information).

There are two Protocols included in the kit: Protocol I and II. The purity and yield of the His tagged protein can often be improved by increasing the wash volume and wash times. Protocol I can be performed easily and adapted to the expected large amount of His tagged protein to be purified. Protocol II includes increasing the wash volume and wash times. Protocol II can be adapted to obtain higher purity and yield of the His tagged protein. We recommend Protocol II when His tagged protein is expected low expression or expressed using mammalian expression systems.

The following protocols are for the isolation of His tagged proteins produced in a 100-mm cell culture dish. The expression level of the His tagged protein may vary. If necessary, adjust the volume of Anti-His tag Beads and Elution Peptide Solution proportionally.

### **Material Preparation**

#### 1. Wash Solution

Dilute Wash Concentrate with 9 times its volume of distilled water.

(e.g. Dilute 0.1 mL of Wash Concentrate with 0.9 mL of distilled water.)

Protocol I; For each Spin Column, prepare 1 mL of Wash Solution.


Protocol II; For each Spin Column, prepare 6 mL of Wash Solution.

#### 2. Elution Peptide Solution

Reconstitute the Elution Peptide with 1 mL of distilled water. If you want to store the reconstituted Elution Peptide, prepare appropriate aliquots (e.g. 45  $\mu$ L x 20 tubes) and store at -20°C. Repeated freezing and thawing is not recommended.

## **Protocol I**

### **Procedure Summary (Purification from mammalian cultured cell lysate)**

### **A. Purification from mammalian cultured cell lysate**

(His tagged protein is not secreted from the cells)

#### **(Lysis of Mammalian Cells)**

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 0.5 mL of Lysis buffer to the cell pellet and vortex.
6. Sonicate the sample for 15 seconds.
7. Incubate the sample for 15 minutes on ice.
8. Remove cell debris by centrifugation at 15,000 x g for 5 minutes at 4°C.

#### **(Purification of His tagged Protein)**

9. Transfer the 0.5 mL of cell lysate (supernatant from step 8) to the Spin Column.
10. Resuspend the Anti-His tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
11. Dispense 20  $\mu$ L Anti-His tag Beads suspension (5  $\mu$ L Beads) into the Spin Column. Screw on the cap.
12. Incubate with gentle end-over-end mixing for 1 hour at 4°C. If the Spin Column does not fit your end-over-end rotator, put it in a suitable tube (e.g. 15 mL centrifuge tube) that fits your end-over-end rotator.
13. Loosen the top cap on the column. Remove the bottom plug. Don't discard the bottom plug. Place the Spin Column in a sampling tube. Centrifuge for 10 seconds. Discard the flow-through (or save for future analysis).
14. Take off the top cap. Keep the bottom plug off. Place the Spin Column in a sampling tube.
15. Add 0.2 mL of Wash Solution to each column. It is not necessary to stir the Spin Column. Centrifuge for 10 seconds. Discard the flow-through. Repeat this step two additional times.

16. Place the Spin Column in a new sampling tube.
  17. For the first elution, screw the bottom plug on tightly. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then screw the top cap on tightly. Tap the tube gently several times. Incubate for 30 minutes at room temperature. Remove the top cap and bottom plug. Centrifuge for 10 seconds.
  18. For the second elution, it is not necessary to place the bottom plugs and top cap on the Spin Column. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then tap the tube gently several times. Incubate for 5 minutes at room temperature. Centrifuge for 10 seconds.
- The two eluates (step 17 and 18) may be pooled in one sampling tube.

## **B. Purification from culture supernatant**


**(His tagged protein is secreted into the culture supernatant)**

1. Collect the culture supernatant from the cell culture dish into a 15 mL centrifuge tube.
2. Centrifuge at 400 x g for 5 minutes to remove cell debris.
3. Transfer the supernatant to a new 15 mL centrifuge tube.
4. Resuspend the Anti-His tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
5. Dispense 20  $\mu$ L Anti-His tag Beads suspension (5  $\mu$ L beads) into the 15 mL centrifuge tube with culture supernatant. Screw on the cap.
6. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
7. Centrifuge at 400 x g for 5 minutes. Discard the supernatant (or save for future analysis), but leave 100-400  $\mu$ L supernatant above the Anti-His tag Beads.
8. Resuspend the Anti-His tag beads in the 100-400  $\mu$ L supernatant by pipetting up and down several times.
9. Transfer the resuspended Anti-His tag Beads in 100-400  $\mu$ L supernatant to the Spin Column.
10. Keep the top cap off. Remove the bottom plug. Don't discard the bottom plug. Place the Spin Column in a sampling tube. Centrifuge for 10 seconds. Discard the flow-through (or save for future analysis). Keep the bottom plug off. Place the Spin Column in a sampling tube.
11. Add 0.2 mL of Wash Solution to each column. It is not necessary to stir the Spin Column. Centrifuge for 10 seconds. Discard the flow-through. Repeat this step two additional times.
12. Place the Spin Column in a new sampling tube.
13. For the first elution, screw the bottom plug on tightly. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then screw the top cap on tightly. Tap the tube gently several times. Incubate for 30 minutes at room temperature. Remove the top cap and bottom plug. Centrifuge for 10 seconds.
14. For the second elution, it is not necessary to place the bottom plugs and top cap on the Spin Column. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then tap the tube gently several times. Incubate for 5 minutes at room temperature. Centrifuge for 10 seconds. The two eluates (step 13 and 14) may be pooled in one sampling tube.

(Note: Steps 11-14 of this protocol are identical to steps 15-18 of the first protocol.)

## Protocol II

### Procedure Summary (Purification from mammalian cultured cell lysate)

### **Purification from mammalian cultured cell lysate**

**(His tagged protein is not secreted from the cells)**

#### **(Lysis of Mammalian Cells)**

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 0.5 mL of Lysis buffer to the cell pellet and vortex.
6. Sonicate the sample for 15 seconds.
7. Incubate the sample for 15 minutes on ice.
8. Remove cell debris by centrifugation at 15,000 x g for 5 minutes at 4°C.

### **(Purification of His tagged Protein)**


9. Transfer the 0.5 mL of cell lysate (supernatant from step 8) to a 1.5 mL sampling tube.
  10. Resuspend the Anti-His tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
  11. Dispense 20  $\mu$ L Anti-His tag Beads suspension (5  $\mu$ L Beads) into the sampling tube (step 9).
  12. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
  13. Centrifuge for 10 seconds. Carefully remove and discard the supernatant (or save for future analysis) using a micropipettor. Don't aspirate beads or disturb bead pellet.
  14. Add 1 mL of Wash Solution to the sampling tube. Inverting the sampling tube several times and centrifuge for 10 seconds. Carefully remove and discard the supernatant using a micropipettor. Don't aspirate beads or disturb bead pellet. Repeat this step four additional times\*.
- \***Note**; Each laboratory is recommended to confirm its own appropriate washing times, because the purity of the His tagged protein may vary depending upon your expression system and so on.
14. Take off the top cap. Keep the bottom plug off. Place the Spin Column in a sampling tube. Don't discard the bottom plug.
  15. Add 0.2 mL of Wash Solution to the sampling tube. Resuspend Anti-His tag Beads by pipetting up and down several times. Transfer the resuspended Anti-His tag Beads to the Spin Column. Place the Spin Column in a new sampling tube.
  16. Centrifuge for 10 seconds. Discard the flow-through.
  17. Place the Spin Column in a new sampling tube.
  18. For the first elution, screw the bottom plug on tightly. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then screw the top cap on tightly. Tap the tube gently several times. Incubate for 30 minutes at room temperature. Remove the top cap and bottom plug. Centrifuge for 10 seconds.
  19. For the second elution, it is not necessary to place the bottom plugs and top cap on the Spin Column. Add 20  $\mu$ L Elution Peptide Solution to the Anti-His tag Beads, then tap the tube gently several times. Incubate for 5 minutes at room temperature. Centrifuge for 10 seconds. The two eluates (step 18 and 19) may be pooled in one sampling tube.

### **Related Products:**

3310	His tagged Protein PURIFICATION KIT	20 purifications
3310A	His tagged Protein PURIFICATION KIT	2 purifications
3311	His tagged Protein PURIFICATION GEL	1 mL
3312	His tagged Protein PURIFICATION GEL	5 mL
3310-205	His tag peptide	2 mg x 5 vials
3305	c-Myc tagged Protein MILD PURIFICATION KIT ver.2	20 purifications
3305A	c-Myc tagged Protein MILD PURIFICATION KIT ver.2	2 purifications
3306	c-Myc tagged Protein MILD PURIFICATION GEL	1 mL
3307	c-Myc tagged Protein MILD PURIFICATION GEL	5 mL
3300-205	c-Myc tag peptide	1 mg x 5 vials
3320	HA tagged Protein PURIFICATION KIT	20 purifications
3320A	HA tagged Protein PURIFICATION KIT	2 purifications
3315	V5 tagged Protein PURIFICATION KIT	20 purifications
3315A	V5 tagged Protein PURIFICATION KIT	2 purifications

## Example of Purification Results 1

### Purification of C-terminus His tagged protein X from culture supernatant

### SDS-PAGE (Coomassie Brilliant Blue Staining)

Stable transfectant of CHO (Chinese Hamster Ovary) cells expressing C-terminus His tagged protein were cultured for 7 days in DMEM medium containing 10% Fetal bovine serum. C-terminus His tagged protein was purified from 0.5 mL of cultured medium according to the preceding protocol I-B. For comparison, elution was carried out not only with peptide solution but also with acid solution and alkaline solution. Each purification was conducted with the same amount of Anti-His tag Beads (5  $\mu$ L) and the same amount of elution solution (20  $\mu$ L x 2 times and then pooled).

Peptide elution solution: neutral pH


0.1 M Glycine-HCl: pH 3.0 (Neutralize the elution immediately with 1 M Tris-HCl, pH 8.0)  
(Acid elution solution)

0.1 M  $\text{NH}_3$ : pH 11.3 (Neutralize the elution immediately with 1 N acetic acid)  
(Alkali elution solution)

## Example of Purification Results 2

### Purification and enzymatic activity of N-terminus His tagged $\beta$ -galactosidase

### SDS-PAGE (Coomassie Brilliant Blue Staining) & X-gal staining

Human embryonic kidney cells (293T) were transfected with pcDNA-His- $\beta$ -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (1 mL/100-mm dish) and purified according to the preceding protocol II. For comparison, elution was carried out not only with peptide solution but also with acid solution and alkaline solution. Each purification was conducted with the same amount of Anti-His tag Beads (5  $\mu$ L) and the same amount of elution solution (20  $\mu$ L x 2 times and then pooled).

Peptide elution solution: neutral pH


0.1 M Glycine-HCl: pH 3.0 (Neutralize the elution immediately with 1 M Tris-HCl, pH 8.0)  
(Acid elution solution)

0.1 M  $\text{NH}_3$ : pH 11.3 (Neutralize the elution immediately with 1 N acetic acid)  
(Alkali elution solution)

Enzymatic activity of each purification was performed using standard X-gal staining method.

## Example of Purification Results 3

### Purification of N- and C-terminus His tagged protein Z from transformed *E. coli*

#### SDS-PAGE (Coomassie Brilliant Blue Staining)

The transformed competent *E. coli* BL21(DE3) cells were induced by the addition of IPTG in the culture medium (1 mL) to express the N- and C-terminus His tagged protein Z. After an induction period, the cells were pelleted by centrifugation and then resuspended in 0.3 mL of the Lysis buffer. His tagged protein Z was purified according to the preceding protocol I-A. For comparison, elution was carried out not only with peptide solution but also with Laemmli SDS-PAGE buffer. Each purification was conducted with the same amount of Anti-His tag Beads (5  $\mu$ L) and the same amount of elution solution (20  $\mu$ L x 2 times and then pooled).

## Additional Information

Several reagents were examined whether or not they were suitable for use with the His tagged Protein PURIFICATION KIT. For example, RIPA buffer could be used for preparation of cell lysate. The results are listed below.

### Chaotropic agents

	<b>Urea</b>	<b>1 M</b>	<b>Yes</b>
	<b>Guanidine-HCl</b>	<b>1 M</b>	<b>No</b>

### Reducing agents

	<b>DTT</b>	<b>10 mM</b>	<b>Yes</b>
	<b>2-Mercaptoethanol</b>	<b>10 mM</b>	<b>Yes</b>

### Surfactants

<b>Nonionic</b>	<b>Tween-20</b>	<b>1%</b>	<b>Yes</b>
	<b>Tween-20</b>	<b>5%</b>	<b>No</b>
	<b>Triton X-100</b>	<b>5%</b>	<b>Yes</b>
	<b>NP-40</b>	<b>1%</b>	<b>Yes</b>
	<b>Digitonin</b>	<b>1%</b>	<b>Yes</b>
	<b>n-Octyl-beta-D-gulcoside</b>	<b>1%</b>	<b>Yes</b>
<b>Zwitterionic</b>	<b>CHAPS</b>	<b>1%</b>	<b>Yes</b>
	<b>CHAPSO</b>	<b>1%</b>	<b>Yes</b>
<b>Anionic</b>	<b>SDS</b>	<b>0.1%</b>	<b>Yes</b>
	<b>Sodium Deoxycholate</b>	<b>0.5%</b>	<b>Yes</b>

### Others

	<b>NaCl</b>	<b>1 M</b>	<b>Yes</b>
	<b>Glycerol</b>	<b>10%</b>	<b>Yes</b>
	<b>EDTA</b>	<b>10 mM</b>	<b>Yes</b>

The “Yes” indicates the reagents can be used in the Lysis buffer for this kit up to the indicated concentration. The “No” indicates the reagents can not be used in the Lysis buffer for this kit at the indicated concentration.

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