CODE No. 4810

MBL

Printed: June 15, 2004
Revised: September 1, 2011

For Research Use Only. Not for use in diagnostic procedures.

Detection of caspase activity

APOPCYTO

Caspase-9 Colorimetric Assay Kit

CODE No. 4810
100 assays

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
URL <https://ruo.mbl.co.jp>
e-mail support@mbi.co.jp, TEL 052-238-1904

Before use, thoroughly read these Instructions.

Description

Caspase is a member of the cysteine aspartic acid-specific protease family, which is activated by a variety of signals of death receptor ligation, DNA damages, serum starvation and stresses *etc.* Active caspase recognizes a lot of several molecules as substrates to cleave them, occurring to biological events corresponding to the apoptosis. For example, ICAD (inhibitor of caspase-activated deoxyribonuclease) is inactivated and CAD (caspase-activated deoxyribonuclease) is indirectly activated by caspase-3, and it is related to chromatin fragmentation for nucleosome units. Caspase recognizes several structural proteins as a substrate to cleave them, and the cleavage is associated with the unique apoptosis cell morphology of chromatin condensation, nucleus fragmentation and cytoplasmic integrity. Active caspase could act in apoptosis process, so that detection of caspase activity is necessary for apoptosis research frequently. Though there are numerous ways to measure of caspase activity, colorimetric or fluorometric substrates are used popularly.

Assay Principle

Active caspase recognizes 4 amino acid residues in a substrate molecule, which are from the aspartate residue next to N-terminal 3 amino acid residues, and specifically cleaves at the C-terminal side of the aspartate residue. So this characterization is useful for detecting caspase activity by using synthetic substrates. In this method, 4 amino acid sequences is labeled with *p*NA (*p*-nitroanilide), MCA (4-Methyl Coumaryl 7-amide) or AFC (7-amino-4-trifluoromethyl coumarin) at the C-terminal side. Free *p*NA, AMC(7-amino-4-methyl-coumarin) or AFC is released from the labeled synthetic substrate on cleavage by active caspase, and monitored by a 96 well microplate reader (*p*NA: absorbance wavelength 400 or 405 nm, AMC: excitation wavelength 380 nm, fluorescence wavelength 460 nm, AFC: excitation wavelength 400 nm, fluorescence wavelength 505 nm). For example, if caspase is active in sample, tetrapeptide-*p*NA as a synthetic substrate is cleaved and the wavelength is shifted as shown in Fig.1. Because the amount of free *p*NA is proportional to the amount of caspase activity present in the sample, caspase activity can be calculated by monitoring of the optical density of free *p*NA at wavelength 405nm.

APOPCYTO Caspase-9 Colorimetric Assay Kit is provided to detect caspase-9 activity in cell extract with LEHD-*p*NA as a substrate, which LEHD sequence is recognized by active caspase-9 selectively. APOPCYTO Caspase-9 Colorimetric Assay Kit is composed of reagents necessary for the detection, which are Cell Lysis buffer, Reaction buffer and so on. Caspase-9 activity can be detected easily and quickly using APOPCYTO Caspase-9 Colorimetric Assay Kit.



Fig. 1 Comparison of UV-Absorption Spectrum of Peptidyl-*p*NA with that of *p*NA

Intended Use

For research use only. Not for use in *in vitro* diagnostic procedures for clinical diagnosis.

Product Components

Materials	Quantity	Color
Cell Lysis Buffer	50 mL × 1 vial	Clear
2X Reaction Buffer	2 mL × 3 vials	Blue
1 M DTT	500 μL × 1 vial	Yellow
Substrate LEHD- <i>p</i> NA (10 mM)	550 μL × 1 vial	Brown
Inhibitor LEHD-FMK (1 mM)	55 μL × 1 vial	Orange
100 mM <i>p</i> NA (<i>p</i> -nitroanilide)	200 μL × 1 vial	Green

Expiration

Please see the label of this kit.

Storage Conditions

Store at -20°C. After thaw, the Cell Lysis Buffer and 2X Reaction Buffer can be stored at 4°C.

Material to Be Supplied by the User

- microcentrifuge (For harvest cells)
- 96-well microplate reader
(For measure the absorbance in the well at 400-405 nm)
- micropipette
- 96-well microplate (flat-bottom, clear polystyrene)
- microcentrifuge tube

Note

- When caspase-9 activity is not measured immediately after cell preparations, store the cell lysates or cell pellets at -20°C and use them within 1 week.
- Store caspase-9 substrate away from light.
- Ensure that DTT is added to 2X Reaction Buffer to obtain a final concentration of 10 mM before use.

Preparation of reagents

- 1) Add DTT to 2X Reaction Buffer to obtain final concentration of 10 mM just before use. (For example, 10 μ L of 1 M DTT to 990 μ L of 2X Reaction Buffer)
- 2) Preparation of *p*NA (*p*-nitroanilide) standards
 - i) Dilute 100 mM *p*NA (*p*-nitroanilide) to 5 mM *p*NA in Cell Lysis Buffer (5 μ L of 100 mM *p*NA is added to 95 μ L of Cell Lysis Buffer).
 - ii) Dilute 5 mM *p*NA in Cell Lysis Buffer according to the following table to give these final concentrations.

<i>p</i>NA concentration	<i>p</i>NA	Cell Lysis Buffer
500 μ M (50 nmole)	50 μ L of 5 mM <i>p</i> NA	450 μ L
250 μ M (25 nmole)	250 μ L of 500 μ M <i>p</i> NA	250 μ L
125 μ M (12.5 nmole)	250 μ L of 250 μ M <i>p</i> NA	250 μ L
62.5 μ M (6.25 nmole)	250 μ L of 125 μ M <i>p</i> NA	250 μ L
31.25 μ M (3.125 nmole)	250 μ L of 62.5 μ M <i>p</i> NA	250 μ L
15.625 μ M (1.5625 nmole)	250 μ L of 31.25 μ M <i>p</i> NA	250 μ L
0 μ M (0 nmole)		250 μ L

Caspase Assay Protocol

1. Induce apoptosis in cells using the desired method, Concurrently incubate a control culture without induction.
2. Count the cells and harvest the $1-5 \times 10^6$ cells by centrifugation at $400 \times g$ for 5 minutes.
3. Resuspend the cells in ice-cold Cell Lysis Buffer (50-500 μ L), and incubate the cells on ice for 10 minutes.
4. Centrifuge the cell lysates at $10,000 \times g$ for 5 minutes at 4°C to precipitate cellular debris. Transfer the supernatants (cell extracts) to new microcentrifuge tubes and put on ice.
If necessary, total protein concentration of the cell extracts may be measured by standard protein assay methods, and adjust protein concentration to 100-200 μg / 50 μL in Cell Lysis Buffer.
5. Add 50 μL of 2X Reaction Buffer containing 10 mM DTT to each well of a 96-well microplate.
10 mM DTT should be added to 2X Reaction Buffer just before use.
6. Add 50 μL of cell lysates to each well. Prepare wells which were added 50 μL of Cell Lysis Buffer in stead of cell lysates to measure blank absorbance.
7. [Optional method]
To verify that the signal detected by the kit is due to protease activity, incubate an induced sample with LEHD-FMK before adding substrate.
8. Add 50 μL of cell extracts to 50 μL of 2X Reaction Buffer (containing 10 mM DTT) and 1 μL of 1 mM LEHD-FMK, and sequentially add 5 μL of Caspase-9 Substrate.
9. Cover the plate with a plate sealer and incubate at 37°C for 2-4 hours (or up to overnight maximum).
10. Add 100 μL of *p*NA (*p*-nitroanilide) standards to empty wells.
11. Measure the absorbance in the wells at 400 or 405 nm. Calculate the specific activity (SA) of caspase-9 present in each sample using following formula:
 - i) Construct a standard curve using the absorbance of *p*NA (*p*-nitroanilide) standards.
 - ii) Calculate the relative absorbance (*A_r*) as:
 $A_a = (\text{induced apoptosis sample } A_{405}) - (\text{blank } A_{405})$
 $A_n = (\text{negative sample } A_{405}) - (\text{blank } A_{405})$
 $A_r = A_a - A_n$
 - iii) Calculate liberated free *p*NA concentration (*B* μM) equal to *A_r* by *p*NA standard curve.
 - iv) Calculate caspase-9 activity (*C*) present in each sample as:
 $C = \text{nmole } p\text{NA liberated per hour}$
 $= B \times 0.1 \text{ mL} / \text{incubation hour}$
 $= 0.1B / \text{incubation hour}$

- v) If total protein concentration of the cell extracts is D mg/mL, Calculate the specific activity (SA) of caspase-9 present in each sample as:

$$SA = C \text{ per } \mu\text{g protein}$$

$$= (0.1B / \text{incubation hour}) / (D \text{ mg/mL} \times 0.1 \text{ mL})$$

$$= B / (D \cdot \text{incubation hour})$$

Assay Example


Fig.2 Caspase-9 activity in CH-11 treated Jurkat

After Jurkat cells were incubated with 100 ng/mL of anti-Fas monoclonal antibody (CH-11) in the presence or absence of caspase-9 inhibitor LEHD-FMK at 37°C for 4 hours, caspase-9 activity was measured by APOPCYTO Caspase-9 Colorimetric Assay Kit.

References

- 1) Wolf, B. B., *et al.*, *J. Biol. Chem.* **274**, 20049-20052 (1999)
- 2) Thornberry, N. A., *et al.*, *J. Biol. Chem.* **272**, 17907-17911 (1997)
- 3) Walker, N. P., *et al.*, *Cell* **78**, 343-352 (1994)
- 4) Sleath, P. R., *et al.*, *J. Biol. Chem.* **265**, 14526-14528 (1990)

Manufacturer

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL <https://ruo.mbl.co.jp>

e-mail support@mbi.co.jp, TEL 052-238-1904

Before use, thoroughly read these Instructions.

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