CODE No. 7620

MBL

Printed: February 10, 1999
Revised: July 19, 1999
Revised: May 1, 2001
Revised: April 2, 2004
Revised: March 20, 2006
Revised: August 01, 2011

For Research Use Only. Not for use in diagnostic procedures.

Quantitative test kit for Human IL-18

Human IL-18 ELISA Kit

CODE No. 7620 96 wells

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
URL <https://ruo.mbl.co.jp>
e-mail support@mbi.co.jp, TEL 052-238-1904

CONTENTS

1. Intended Use -----	1
2. Summary and Explanation-----	1
3. Principle -----	2
4. Materials provided -----	2
5. Storage and Stability -----	2
6. Species and cross reactivity-----	3
7. Materials and equipment required -----	3
8. Precautions -----	3
9. Procedure -----	4
10. Performance Characteristics -----	8
11. Related products-----	10
12. References -----	11

CODE No. 7620

Before use, thoroughly read these Instructions.

Intended Use

The **Human IL-18 ELISA Kit** is based on sandwich ELISA and is capable of measuring human IL-18.

For Research Use Only. Not for use in diagnostic procedures.

Summary and Explanation

Interleukin 18 (IL-18) is an 18 kDa cytokine which is identified as a costimulatory factor for production of interferon- γ (IFN- γ) in response to toxic shock. It shares functional similarities with IL-12. IL-18 is synthesized as a precursor 24 kDa molecule without a signal peptide and must be cleaved to produce an active molecule. IL-1 β converting enzyme (ICE, caspase-1) cleaves pro-IL-18 at aspartic acid in the P1 position, producing the mature, bioactive peptide that is readily released from the cells. It has been reported that IL-18 is produced from dendritic cells, activated macrophages, Kupffer cells, keratinocytes, intestinal epithelial cells, osteoblasts, adrenal cortex cells and murine diencephalons.

IFN- γ is produced by activated T and NK cells and plays critical roles in the defense against microbial pathogens. IFN- γ activates macrophages, enhances NK activity and B cell maturation, proliferation and Ig secretion, induces MHC class I and II antigens expression, and inhibits osteoclast activation.

IL-18 acts on T helper 1-type (Th1) cells, and in combination with IL-12 strongly induces production of IFN- γ by these cells. Pleiotropic effects of IL-18 have also been reported, including enhancement production of IFN- γ and GM-CSF in peripheral blood mononuclear cells, production of Th1 cytokines, IL-2, GM-CSF and IFN- γ in T cells, enhancement of Fas ligand expression by Th1 cells.

The "Human IL-18 ELISA Kit" is a useful reagent for specifically measuring human IL-18 with high sensitivity by ELISA.

This Kit does not detect 1 ng/mL of various cytokines, such as human IFN- α , IFN- γ , IL-1 β , IL-4, IL-5, IL-6, IL-10, IL-12, GM-CSF and murine IL-18. The results were all below the detection limit of 12.5 pg/mL.

Principle

The **Human IL-18 ELISA Kit** measures human IL-18 by sandwich ELISA. This assay uses two monoclonal antibodies against two different epitopes of human IL-18. Samples to be measured or standards are incubated in the microwells coated with anti-human IL-18 monoclonal antibody. After washing, the peroxidase conjugated anti-human IL-18 monoclonal antibody is added into the microwells and incubated. After another washing, the substrate reagent is mixed with the chromogen and allowed to incubate for an additional period of time. An acid solution is added to each microwell to terminate the enzyme reaction and to stabilize the color development. The optical density (O.D.) of each microwell is measured at 450 nm using a microplate reader.

The concentration of human IL-18 is calibrated from a dose response curve based on the reference standards.

Materials provided

Each kit contains:

Name	Materials	Quantity (for 1 kit)
Microwells	microwell strips coated with anti-human IL-18 antibody	8-well strip × 12 strips
Human IL-18 calibrator	human IL-18 (lyophilized powder)	2 vials
Conjugate reagent	Peroxidase conjugated anti-human IL-18 antibody (× 101)	0.2 mL × 1 vial
Conjugate diluent	Buffer for diluting Conjugate reagent (ready-for-use)	24 mL × 1 bottle
Assay diluent	Buffer for diluting samples* (ready-for-use)	30 mL × 1 bottle
Wash concentrate	Buffer for washing microwells (× 10)	100 mL × 1 bottle
Substrate reagent	TMB/H ₂ O ₂ solution (ready-for-use)	20 mL × 1 bottle
Stop solution	0.5 M H ₂ SO ₄ (ready-for-use)	20 mL × 1 bottle

* Use the **Assay diluent** as 0 pg/mL of IL-18 standard.

Storage and Stability

All kit components must be stored at 2-8°C. All reagents are stable for 12 months after manufacturing when stored at the indicated conditions.

Species cross reactivity

Species	Human	Mouse	Rat	Monkey
Samples	serum, plasma, cell culture supernatant, urine*	recombinant	Not Tested	Not Tested**
Reactivity	+	-		

* This kit was used in reference 5).

** This kit was used in reference 9).

Materials and equipment required

- Microplate reader (450 nm)
- Automatic plate washer or wash bottle
- Adjustable micropipette
- Multichannel micropipette
- 96-well polyvinyl plate
- Microplate holder
- Uncoated microwell strips (for using an automatic plate washer)
- Disposable reagent vessels
- Paper towels
- Distilled water

Precautions

1. Allow all the components to come to room temperature (20-25°C) before use.
2. **Microwells** that are not immediately required should be returned to the ziplock pouch, which must be carefully resealed to avoid moisture absorption.
3. Fresh samples should be used. Aliquot each sample and store below -20°C if necessary. Avoid repeated freezing and thawing. Never store the samples at 4°C, as samples might be affected by storage of this temperature.
4. When **Wash concentrate** is stored at 2-8°C, some precipitation or turbidity may appear. However it does not affect the reagent efficiency. Allow the **Wash concentrate** to come to room temperature and mix to re-dissolve the precipitate.
5. **Assay diluent** contains sodium azide (0.09%) as preservative. Azide may react with copper or lead in plumbing system to form explosive metal azide. Therefore, always flush with plenty of water into a drain when disposing materials containing azide.
6. **Stop Solution** is 0.5 M sulfuric acid. As it is a corrosive material, protect eyes and skin and handle with care.
7. This kit is intended for research use only. Not for use in diagnostic procedures.

Procedure

◆ Preparation of Reagents

1. Wash solution

Prepare “**Wash solution**” by diluting **Wash concentrate, 1:10** with distilled water prior to use. (e.g., add 100 mL of Wash Concentrate to 900 mL of distilled water.) This “**Wash solution**” is stable for 2 weeks at 4°C.

2. Conjugate solution

Prepare “**Conjugate solution**” by diluting the **Conjugate reagent, 1:101** with the **Conjugate diluent** prior to use. (e.g., add 10 µL of Conjugate reagent to 1,000 µL of Conjugate diluent.)

* Prepare only a sufficient amount of the “**Conjugate solution**” for the assay because the diluted Conjugate reagent is not stable.

* Use disposable new pipette and vessel to avoid contamination of microbe.

3. Standards

Reconstitute “**Standards**” by dissolving lyophilized **Human IL-18 calibrator** with the **Assay diluent**. The volume of **Assay diluent**, and precise procedures are indicated in the attached lot specific document “*Preparation of Standards*”.

* If storage is needed after reconstitution, prepare appropriate aliquots and freeze them below -20°C. Avoid repeated freezing and thawing.

4. Other reagents are ready for use.

◆ Preparation of samples

1. Sampling

Each laboratory is recommended to establish its own sampling procedure. The following sampling procedures have been reported.

1) Serum

Blood for measurement of IL-18 was collected in Venoject®-II Tube AUTOSEP tube (TERUMO, Japan; cat# VP-AS109K). The blood was allowed to clot at room temperature for exactly 60 ± 10 minutes. Thereafter, the serum was separated by centrifugation twice for 10 minutes at 3,000 rpm at 4°C, and then transferred to a new plastic test tube.

CODE No. 7620

2) EDTA Plasma

Blood for measurement of IL-18 was collected in Venoject[®]-II Tube (TERUMO, Japan; cat# VP-DK050K). Within 30 minutes of blood collection, the plasma was separated by centrifugation twice for 10 minutes at 3,000 rpm at 4°C, and then transferred to a new plastic test tube.

2. Dilution

1) Dilute samples

Dilute each sample **1:5** with **Assay diluent**.

Example: human serum or EDTA plasma

Adding 50 µL of sample to 200 µL of **Assay diluent**.

2) Add 150 µL of prepared samples and **Standards** to 96-well polyvinyl plate as the same order of assay run.

3. Storage

Fresh samples should be used. Aliquot each sample into new plastic tube, and store below -20°C if necessary. Avoid repeated freezing and thawing.

◆ Assay procedure

Duplicate assay is recommended. A standard curve must be run with each assay.

STEP 1. (Sample incubation)

1) Add 150 µL of prepared samples and **Standard** to a transient 96-well polyvinyl plate as the same order of actual assay run. Then, transfer 100 µL of each sample to **Microwells** simultaneously using multichannel pipette.

Because the reaction starts on the transfer to the well, this transfer should be completed as quickly as possible.

2) Cover the plate and incubate for 60 minutes at room temperature (20-25°C).

STEP 2. (Washing)

Manual wash: Aspirate or discard the well contents, then fill each well with **Wash solution** using wash bottle. Aspirate or discard the **Wash solution** in the wells. Repeat this step another 3 times. Tap the plate on a paper towel to remove any remaining **Wash solution**.

Autowasher: Wash 4 times with **Wash solution**. Tap the plate on a paper towel to remove any remaining **Wash solution**.

* Each laboratory is recommended to confirm its own appropriate washing condition.

* **Wash solution** should be used at room temperature.

STEP 3. (Conjugate incubation)

- 1) Pour **Conjugate solution** into the vessel. After removing wash solution completely, pipette 100 μL of **Conjugate solution** to each well with a multichannel pipette.
- 2) Cover the plate and incubate for 60 minutes at room temperature.

STEP 4. (Washing)

Wash the wells again following the **STEP 2.** procedure.

STEP 5. (Substrate incubation)

- 1) After the washing, pour **Substrate reagent** (ready-for-use) into vessel. Add 100 μL of **Substrate reagent** to each well.

Because the enzyme reaction starts on the transfer to the well, this transfer should be completed as quickly as possible. For testing many samples, prepare the **Substrate** in a reagent vessel, then pour it into each well using multichannel pipette. This vessel should be different from the one that was used for pouring **Conjugate**.

* **Substrate reagent** should be used at room temperature.

* Use disposable new pipette and vessel, as **Substrate reagent** is easily oxidized by metal ions and can be contaminated by microbes.

* Once **Substrate reagent** is poured into the vessel, do not return to the bottle.

- 2) Cover the plate and incubate for 30 minutes at room temperature.

STEP 6. (Stopping reaction)

Pour **Stop solution** (ready-for-use) into the vessel. Add 100 μL of **Stop solution** to each well with a multichannel pipette.

◆ **Reading**

Read the absorbance of each well at 450 nm using a plate reader. If a dual wavelength plate reader is available, set the test wavelength at 450 nm and the reference at 620 nm.

* Reading should be done within 30 minutes after stopping the reaction.


* Ensure that the back of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading.

◆ **Calculation of results**

- 1) Calculate the mean absorbance value of each **IL-18 standard**.
- 2) Plot on the semi-log graph paper and construct a calibration curve.
[Absorbance on the vertical axis, concentration (pg/mL) on the horizontal axis]
- 3) Read the pg/mL corresponding to the absorbance of sample from the standard curve.
Report the IL-18 concentration of samples by multiplying dilution factor. (e.g. human serum: ×5)

* If absorbance of sample exceeds the one of the 1,000 pg/mL standard, dilute the sample and measure again.

◆ **Example of standard curve**

◆ **Concentration of Human IL-18 in normal human sample**

Serum samples from 46 healthy donors were measured with the Human IL-18 ELISA Kit. The results were as follows. (Calculated as described in “◆ Calculation of results”.)

Maximum	=	257.8 pg/mL
Minimum	=	36.1 pg/mL
Mean	=	126.0 pg/mL
SD	=	44.5 pg/mL
Mean + 3 SD	=	259.4 pg/mL

Performance Characteristics

◆ Sensitivity

The sensitivity of the assay is 12.5 pg/mL.

The minimum detection limit estimated by serial dilution was 12.5 pg/mL since the mean + 2SD of the 6.25 pg/mL was lower than the mean - 2SD of the 12.5 pg/mL.

◆ Reproducibility

1. Intra-assay

Intra-assay reproducibility was determined by assaying the sera 8 times.

IL-18 concentrations of the serum samples were calculated as described in

“**◆ Calculation of results**”.

Sample	serum 1	serum 2	serum 3	serum 4	serum 5
Number of determinations	8	8	8	8	8
Mean (pg/mL)	2765.5	600.7	345.5	136.1	69.7
C.V. (%)	5.03	4.93	10.80	5.61	9.92

2. Inter-assay

Inter-assay reproducibility was determined by 5 independent assays.

IL-18 concentrations of the serum samples were calculated as described in

“**◆ Calculation of results**”.

Sample	serum 1	serum 2	serum 3	serum 4	serum 5
Number of determinations	5	5	5	5	5
Mean (pg/mL)	2621.1	773.2	615.1	236.4	160.1
C.V. (%)	10.07	8.54	5.21	7.60	6.25

* From six replicates of each serum sample in five separate assays.

◆ **Recovery test**

Recombinant human IL-18 was added to samples at different concentrations.

IL-18 concentrations of the serum samples were calculated as described in

“◆ **Calculation of results**”.

serum 1

(A) additional rhIL-18 (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	940.8	-	-
684.9	1583.0	642.2	94
1346.9	2179.0	1238.2	92
2354.1	3528.7	2587.9	110

serum 2

(A) additional rhIL-18 (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	868.9	-	-
684.9	1418.6	549.8	80
1346.9	2035.9	1167.0	87
2354.1	3431.3	2564.2	109

serum 3


(A) additional rhIL-18 (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	629.3	-	-
684.9	1334.7	705.4	103
1346.9	1913.4	1284.1	95
2354.1	3282.3	2653.1	113

* rhIL-18 is abbreviation of recombinant human IL-18.

◆ Dilution test

Serum samples were diluted with **Assay diluent**.

Each serum was serially diluted from 1/1 to 1/32 with **Assay diluent**, then IL-18 concentrations of the serum samples were calculated as described in “◆ Calculation of results”.

Related products

- 7625 Mouse IL-18 ELISA Kit
- D043-3 anti-Human IL-18 (25-2G)
- D044-3 anti-Human IL-18 (125-2H)
- D045-3 anti-Human IL-18 (159-12B)
- D045-6 Biotin labeled anti-Human IL-18 (159-12B)
- PM014 anti-Human IL-18 (poly)
- D046-3 anti-Mouse IL-18 (39-3F)
- D047-3 anti-Mouse IL-18 (74)
- D048-3 anti-Mouse IL-18 (93-10C)
- D048-6 Biotin labeled anti-Mouse IL-18 (93-10C)
- M156-3 anti-Human pro-IL-18 (43A11)
- M157-3 anti-Rat IL-18 (21A12)

CODE No. 7620

- M158-3 anti-Rat IL-18 (91D8)
- M159-3 anti-Human IL-18 Receptor 1 (44G6)
- M163-3 anti-Mouse IL-18 Receptor 1 (33A11)
- M166-3 anti-Mouse IL-18 Receptor 1 (64G4)
- B001-5 Recombinant Human IL-18
- B003-5 Recombinant Human IL-18 (without BSA)
- B002-5 Recombinant Mouse IL-18
- B004-5 Recombinant Mouse IL-18 (without BSA)

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