For Research Use Only. Not for use in diagnostic procedures.

Detection System for Human sections

Histostar™ (Ms + Rb)

Immunostaining Reagent, Anti-Mouse and Anti-Rabbit

CODE No. 8460A, 8460, 8461

Histostar™ (Rb)

Immunostaining Reagent, Anti-Rabbit

CODE No. 8466A, 8466, 8467


1. INTRODUCTION

Immunohistochemistry (IHC) is an efficient means to understand antigens distribution and localization in cells or tissue sections under the microscope. IHC is widely used in clinical diagnosis of disease and basic research purposes. IHC uses monoclonal or polyclonal antibodies to determine the tissue distribution of an antigen of interest. The site of antibody binding is visualized by a marker such as fluorescent dye, enzyme, radioactive element or colloidal gold that is directly linked to the primary antibody or to an appropriate secondary antibody. In immunoenzyme technique, frequently used streptavidin-biotin complex (SAB) methods or avidin-biotin complex (ABC) method, which need additional two steps containing both Biotinylated secondary antibody reaction and peroxidase labeled streptavidin reaction after first antibody reaction.

Histostar™ is a two-step IHC staining reagent. This reagent is based on an horseradish peroxidase (HRP) labeled polymer which is conjugated with secondary antibodies. It increases the number of available HRP at the antigenic site, thus increasing their reactivity with the chromogen. **Histostar™** is an extremely sensitive reagent and save time. Because **Histostar™** does not contain avidin and biotin, nonspecific staining as a result of endogenous biotin is eliminated.

2. PRINCIPLE

The antigen / antibody / **Histostar™** complex can be prepared by allowing the reagent to react with a primary antibody bound to the antigen on tissue section. The enzymatic activity of this complex results in a colored dependent, thus staining the antigen site.

3. PRESENTATION

Liquid. Ready for use.

Histostar™ is stored in MOPS (3-Morpholinopropanesulfonic acid) buffer (pH 6.5) containing stabilizer and antibacterial agent.

Histostar™ (Ms + Rb) (Immunostaining Reagent, Anti-Mouse and Anti-Rabbit) is the labeled polymer prepared by combining amino acid polymers with peroxidase and goat anti-Mouse Ig and goat anti-Rabbit Ig which are reduced to Fab'. IgG fraction purified from immunized goat serum is digested to prepare F(ab')₂. Antigen-specific F(ab')₂ is affinity-purified with the antigen. Solid-phase absorption is carried out with human serum protein. Peroxidase-labeled amino acid polymer is conjugated to Fab' obtained by reducing F(ab')₂.

Histostar™ (Rb) (Immunostaining Reagent, Anti-Rabbit) is the labeled polymer prepared by combining amino acid polymers with peroxidase and goat anti-Rabbit Ig which are reduced to Fab'. IgG fraction purified from immunized goat serum is digested to prepare F(ab')₂. Antigen-specific F(ab')₂ is affinity-purified with the antigen. Solid-phase absorption is carried out with human serum protein. Peroxidase-labeled amino acid polymer is conjugated to Fab' obtained by reducing F(ab')₂.

Products name	Code no.	Quantity
Histostar™ (Ms + Rb)	8460A	1 mL × 1 vial (10 tests)
	8460	15 mL × 1 vial (150 tests)
	8461	15 mL × 3 vials (450 tests)
Histostar™ (Rb)	8466A	1 mL × 1 vial (10 tests)
	8466	15 mL × 1 vial (150 tests)
	8467	15 mL × 3 vials (450 tests)

4. STORAGE and STABILITY

Histostar™ reagents are to be stored in a dark place at 2-8°C. Do not freeze. The expiration date is printed on the reagent vial.

5. INTENDED USE

For Research Use Only. Not for use in diagnostic procedures.

Histostar™ (Ms + Rb) and **Histostar™** (Rb) are designed to allow Immunohistochemical studies using a user-supplied Mouse primary antibody and Rabbit primary antibody, or Rabbit primary antibody, respectively. **Histostar™** reagents are available for formalin fixed paraffin-embedded human tissue sections.

6. PRECAUTION WHILE USING OR HANDLING

1. Before using this reagent, please read these instructions.
2. Do not use reagents after the expiration date.
3. Specimens, before and after fixation, and all other materials exposed to them, should be handled like biohazardous materials with proper precautions.
4. Inhalation or ingestion of the highly allergic fixative formaldehyde is harmful. Wear protective mask. If swallowed, induce vomiting. If skin or eye contact occurs, wash thoroughly with water.
5. Organic reagents are flammable. Do not use near open flame.
6. Never pipette reagents by mouth and avoid their contact with skin, mucous membranes and clothes.
7. Avoid microbial contamination of reagents as incorrect result may occur.
8. Avoid splashing of reagents or generation of aerosols.
9. For research use only. Not for diagnostic use.

7. STAINING PROCEDURES

■ Reagents and Materials required but not provided

- Primary antibody
- Xylene
- 100% ethanol
- 95% ethanol
- Counter staining solution
- 3% solution of hydrogen peroxide in PBS or 0.3% hydrogen peroxide in methanol
- Phosphate buffered saline (PBS)
- Adhesive for tissue section (0.02% poly-L-Lysine, silane or the like)
- Chromogen/substrate reagent
Histostar™ DAB substrate solution (MBL; Code no. 8469) is recommended.
- Negative control reagent
- Cover slips
- Light microscope
- Distilled water
- Humidified chamber for slide incubation
- Mounting media
- Timer
- Staining racks or Coplin Jars
- Absorbent wipes

■ Specimen preparation

[Paraffin-embedded tissue sections]

Specimens may undergo histological disintegration or antigenic denaturation when subjected to highly concentrated fixative or prolonged fixing. Thus, in order to obtain an optimal fixing, maintaining tissue morphology and antigen activity, tissues which are as fresh as possible and small in size (about 1 cm x 1 cm x 0.5cm) should be used. The fixatives as shown are recommended.

Fixing reagent	Fixing time
10% formalin or buffered formalin	24 – 48 hours
20% formalin	12 – 24 hours

[Frozen tissue sections]

Specimens are embedded in compounds (like OCT compound) and snap-frozen in n-hexane cooled in dry ice-acetone or liquid nitrogen.

■ Section preparation

[Paraffin embedded tissue sections]

The cut sections should be 3-6 µm and placed on slides. When further treatments are to be done such as Antigen Recovery, Heat-Induced Epitope Retrieval (HIER) or trypsin treatment, the glass slides should be coated with an adhesive like 0.02% poly-L-lysine or silane for tissue sections.

[Frozen tissue sections]

The Cryostat sections are mounted on an adhesive (like 0.02% poly-L-Lysine)-coated slides and air-dried well. The sections are fixed in 100% acetone for 10 minutes at room temperature or 4% paraformaldehyde-PBS solution for 10 minutes at 4°C and then stained.

[Control slides]

A positive control slide, negative control slide and reagent control slide are needed and processed in the same way as the unknown specimen slide to interpret staining results.

Protocol

■ Deparaffinization and Rehydration

1. Immerse the slides in Xylene 3 times for 3-5 minutes each.
2. Immerse the slides in 100% Ethanol 2 times for 3-5 minutes each.
3. Immerse the slides in 95% Ethanol 2 times for 3-5 minutes each.
4. Immerse the slides in PBS 3 times for 3-5 minutes each.

■ Recommended Staining Procedures

1. Wipe areas around the sections on the slides carefully to remove excess solution.
2. Immerse each section in 3% H₂O₂ for 10 minutes or 0.3% H₂O₂ in Methanol for 30 minutes at room temperature to block endogenous peroxidase activity.
3. Wash the slides 3 times in PBS for 5 minutes each.
4. Remove the slides from PBS, wipe gently around each section and cover tissues with Blocking reagent (e.g. 20 mM HEPES, pH 7.4, 1% BSA, 135 mM NaCl) for 10-15 minutes to block non-specific staining. Do not wash.
5. Tip off the blocking buffer, wipe gently around each section and cover tissues with optimally diluted primary antibody or negative control reagent. Incubate them sections at room temperature or 4°C. (Follow the instructions for dilute concentration and incubation time designated in the package insert of primary antibody. The optimal concentration and incubation time of primary antibodies should be determined by the investigation. In some cases, further dilution of primary antibodies may be required to prevent overstaining.)
6. Wash the slides 3 times in PBS for 5 minutes each.
7. Wipe gently around each section and cover tissues with **HistostarTM** reagent. Incubate for 30 minutes at room temperature.
8. Wash the slides 3 times in PBS for 5 minutes each.
9. Visualize by reacting for 5-20 minutes with HistostarTM DAB substrate solution. *DAB is a suspect carcinogen and must be handled with care. Always wear gloves.
10. Wash the slides in water for 5 minutes.
11. Counter stain in hematoxylin for 1 minute, wash the slides 3 times in water for 5 minutes each, and then immerse the slides in PBS for 5 minutes. Dehydrate by immersing in Ethanol 3 times for 3 minutes each, followed by immersing in Xylene 3 times for 3 minutes each.
12. Now ready for mounting.

■ Interpretation of results

- Microscopic observation

The slides are examined under a light microscope for a positive reaction. It is necessary to make comparison with three types of the control slides for interpreting staining results.

- Positive control slide

A specimen containing the target antigen which is processed in the same way as the unknown specimen.

- Negative control slide

A specimen not containing the target antigen which is processed in the same way as the unknown specimen.

- Reagent control slide

The control specimen is used and processed in the same way as the test specimen except that negative control reagent is used instead of primary antibody. If the slide is stained, it is probably due to non-specific reaction by non-specific protein binding.

8. LIMITATION

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, or sectioning may produce artifacts or false-negative results.

Results will not be optimal if old or unbuffered fixatives are used, or excessively heated during embedding or during attachment of sections to slides.

False-Positive results may be seen due to nonspecific binding of proteins. Although **Histostar™** does not require the use of blocking reagent separately, in some cases the application of blocking reagent containing an irrelevant protein, prior to incubation with the primary antibody, may be useful for reducing the background.

9. CONDITION FOR USE

Histostar™ is designed for research use only and is not intended for therapeutic or diagnostic purposes. MEDICAL & BIOLOGICAL LABORATORIES sales agents and distributors will take no responsibility for **Histostar™** when used in a way which directly or indirectly violates local regulations or patents. Neither MEDICAL & BIOLOGICAL LABORATORIES nor its sales agents can be held responsible for any patent infringement which may occur as a result of improper use of the product.

10. TROUBLE SHOOTING

Problem	Possible cause	Solution
<input type="checkbox"/> No staining or only weak staining results on the positive control slide and the unknown specimen slide.	<ol style="list-style-type: none"> 1. Drying-out of specimens during staining prior to addition of the reagents. 2. The embedding agent is not suitable, or paraffin is not thoroughly removed from paraffin-embedded sections. 3. Any trace amount of sodium azide present in the buffer inactivates the peroxidase, making the staining impossible. 4. Inadequate incubation of the enzyme and antibody. 	<ol style="list-style-type: none"> 1. Never allow the tissue to dry out. 2. Select a suitable embedding agent or remove paraffin thoroughly from sections embedded 2. Change xylene or ethanol as the case may be. 3. Use sodium azide free buffer solution. 3. Change buffer solution. 4. Change stale chromogen/substrate reagent. 4. Blot off excess solution thoroughly at each stage. 4. Provide sufficient time for reaction with antibody. In particular, primary antibody should be incubated for the time period specified in the insert.
<input type="checkbox"/> The unknown specimen slide is not stained while the positive control slide is stained.	<ol style="list-style-type: none"> 1. Antigen is denatured or masked during fixing or embedding process. 2. Antigen is decomposed by autolysis. 3. Less antigen is present in the sections. 	<ol style="list-style-type: none"> 1. Some antigens are sensitive to fixation or embedding. So use less potent fixative and decrease the fixing time. 1. The pretreatment is required for some tissues, in order to reveal the antigen, such as Antigen Recovery, Heat-Induced Epitope Retrieval (HIER) or trypsin treatment. 2. Use tissues obtained by biopsy or surgery, whenever possible. 3. Prolong the incubation time.
<input type="checkbox"/> The backgrounds are intensively stained in all the slides.	<ol style="list-style-type: none"> 1. Endogenous enzyme activity was not completely blocked. 2. Non-specific binding compositions are found. 3. Autolysis results in excessive antigens isolated in histological solutions. 4. Insufficient removal of paraffin. 5. Insufficient washing of antibody. 6. A high room temperature accelerates enzyme reactions. 7. Drying-out of specimens during staining after of the reagents. 	<ol style="list-style-type: none"> 1. Ensure that the procedure for quenching of endogenous peroxidase is right. 2. Before adding primary antibody, treat with 10% normal goat serum. 3. Obtain fresh tissues whenever available. 4. Change xylene or ethanol as the case may be. 5. Ensure thorough washing of antibody. 6. Keep room temperature at 15 to 25°C. 6. Shorten reaction time. 7. Never allow the tissue to dry out.
<input type="checkbox"/> During the reaction, tissue sections come off from the slides.	<ol style="list-style-type: none"> 1. Some antigens require heat induced antigen retrieval procedure or prolonged reaction time with primary antibody, which may render the sections easily come off. 	<ol style="list-style-type: none"> 1. Mount tissue sections on slides coated with an adhesive such as 0.02% poly-L-lysine or silane.

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