For Research Use Only. Not for use in diagnostic procedures.
--

PURIFICATION that Maintains Protein Activity

c-Myc tagged Protein
MILD PURIFICATION KIT
(MoAb. clone 1G4) ver. 2

CODE No. 3305

Purification with neutral pH (to maintain Activity and Conformation)

Using Spin Column (Rapid and Easy)

High Quality and Efficiency

High-Yield Purification

Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more.

The method of purification is one of the important keys for maintaining protein structure and function. The c-Myc tagged Protein MILD PURIFICATION KIT is designed for the isolation of c-Myc tagged protein from cell lysates and culture supernatant. Severe conditions such as acidic or alkaline elution can denature protein structure while maintaining a neutral pH can preserve protein activity and conformation. MBL has developed the Anti-c-Myc tag Gel to quickly and efficiently purify c-Myc-tagged proteins at neutral pH to maintain protein activity and native conformation. The c-Myc tag peptide competitively elutes the c-Myc tagged protein from Anti-c-Myc tag Beads while minimally interfering with protein function.

The affinity of the anti-c-Myc tag antibody is also very important. Very high affinity antibody cannot be dissociated from the c-Myc tagged protein, and low affinity antibody is not sufficient for binding. MBL has optimized the affinity of the anti-c-Myc tag antibody to permit high yields and efficient purification of functional c-Myc-tagged proteins for biochemical studies. This kit is simple and fast, because all procedures of this kit have been optimized by using a Spin Column.


In ver.2, MBL improved the purification capacity of the Anti-c-Myc tag Beads by using the new monoclonal antibody (Clone 1G4). It is easily possible for scale-up by using a c-Myc tagged Protein MILD PURIFICATION GEL (code 3306, 3307).

Kit Components

Components sufficient for conducting 20 purifications of c-Myc tagged protein from 100-mm culture dishes.

- 1. Anti-c-Myc tag Beads** 25% slurry: 100 μ L beads in 400 μ L total volume in PBS with 0.1% ProClin 150 as preservative
- 2. Elution Peptide** c-Myc tag peptide (EQKLISEEDL), 1 mg in 1 mL PBS after reconstitution
- 3. Spin Columns Sets** 20 columns with pre-inserted bottom plugs and top caps
- 4. Wash Concentrate** 10x concentrate, 6 mL

Procedure Summary

Product Capacity

The purification capacity of the Anti-c-Myc tag Beads varies depending upon the c-Myc tagged protein. For examples, 5 μ L of Anti-c-Myc tag Beads (20 μ L slurry) bound 5 μ g of a c-Myc tagged protein (27 kDa) and eluted 3 μ g of purified protein.

Materials Required but not Provided

1. Microcentrifuge capable of 15,000 x g
2. Sampling tube (1.5 mL or 2 mL)
3. End-over-end rotator
4. PBS or TBS
5. Lysis buffer

Suitable Lysis buffer varies with cell type.

Note: Do not use SDS based RIPA buffer (see Additional Information).

A. Homemade Lysis buffer

20-50 mM Tris-HCl, pH 7.5 or HEPES-KOH, pH 7.5

50-250 mM NaCl

5 mM EDTA

1% NP-40 or Triton X-100

if necessary add Protease Inhibitor Cocktail

(e.g. SIGMA code P8340, PIERCE code 78415).

B. Commercial Reagent

SIGMA CelLytic M Cell Lysis Reagent, code no. C2978

PIERCE M-PER Mammalian Protein Extraction Reagent, code no.78503

Storage

Store for up to 1 year from date of receipt at 4°C. Do not freeze.

Protocols

The following protocol is for the isolation of c-Myc tagged proteins produced in a 100 mm cell culture dish. The expression level of the c-Myc tagged protein may vary. If necessary, adjust the volume of Anti-c-Myc tag Beads and Elution Peptide Solution proportionally.

Material Preparation

1. Wash Solution

Dilute Wash Concentrate with 9 times its volume of distilled water. (e.g. Dilute 0.1 mL of Wash Concentrate with 0.9 mL of distilled water.) For each Spin Column, prepare 1 mL of Wash Solution.

2. Elution Peptide Solution

Reconstitute the Elution Peptide with 1 mL of distilled water. If you want to store the reconstituted Elution Peptide, prepare appropriate aliquots (e.g. 45 μ L x 20 tubes) and store at -20°C. Repeated freezing and thawing is not recommended.

Purification from cell lysate (c-Myc tagged protein is not secreted from the cells)

(Lysis of Mammalian Cells)

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS or TBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 500 μ L of Lysis buffer to the cell pellet and vortex.
6. Incubate the sample for 30 minutes on ice.
7. Remove cell debris by centrifugation at 15,000 x g for 10 minutes at 4°C.

(Purification of c-Myc tagged Protein)

8. Transfer the 500 μ L of cell lysate to the Spin Column.
9. Resuspend the Anti-c-Myc tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
10. Dispense 20 μ L Anti-c-Myc tag Beads suspension (5 μ L beads) into the Spin Column. Screw on the cap.
11. Incubate with gentle end-over-end mixing for 1 hour at 4°C. If the Spin Column does not fit your end-over-end rotator, put it in a suitable tube (e.g. 15 mL centrifuge tube) that fits your end-over-end rotator.
12. Loosen the top cap on the column. Remove the bottom plug. Place the Spin Column in a sampling tube. Centrifuge for 10 seconds. Discard the flow-through (or save for future analysis).
13. Take off the top cap. Keep the bottom plug off. Place the Spin Column in a sampling tube.
14. Add 0.2 mL of Wash Solution to each column. It is not necessary to stir the Spin Column. Centrifuge for 10 seconds. Discard the flow-through. Repeat this step two additional times.
15. Place the Spin Column in a new sampling tube.
16. Add 20 μ L Elution Peptide Solution to the Anti-c-Myc tag Beads. It is not necessary to place the bottom plugs and top cap on the Spin Column. Tap the tube gently several times. Incubate for 5 minutes at 4°C. Centrifuge for 10 seconds.
17. Repeat step 16 again. The two eluates may be pooled in one sampling tube.

Purification from culture supernatant (if the c-Myc tagged protein is secreted into the culture supernatant)

1. Collect the culture supernatant from the cell culture dish into the centrifuge tube.
2. Centrifuge at 400 x g for 5 minutes to remove cell debris.
3. Transfer the supernatant to a new 15 mL centrifuge tube.
4. Resuspend the anti-c-Myc tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
5. Dispense 20 μ L anti-c-Myc tag Beads suspension (5 μ L beads) into the 15 mL centrifuge tube with culture supernatant. Screw on the cap.
6. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
7. Centrifuge at 400 x g for 5 minutes. Discard the supernatant (or save for future analysis), but leave 100-400 μ L supernatant above the Anti-c-Myc tag Beads.
8. Resuspend the Anti-c-Myc tag Beads in the 100-400 μ L supernatant by pipetting up and down several times.
9. Transfer the resuspended Anti-c-Myc tag Beads in 100-400 μ L supernatant to the Spin Column.
10. Keep the top cap off. Remove the bottom plug. Place the Spin Column in a sampling tube. Centrifuge for 10 seconds. Discard the flow-through (or save for future analysis).

11. Keep the bottom plug off. Place the Spin Column in a sampling tube.
12. Add 0.2 mL of Wash Solution to each column. It is not necessary to stir the Spin Column. Centrifuge for 10 seconds. Discard the flow-through. Repeat this step two additional times.
13. Place the Spin Column in a new sampling tube.
14. Add 20 μ L Elution Peptide Solution to the Anti-c-Myc tag Beads. It is not necessary to place the bottom plugs and top cap on the Spin Column. Tap the tube gently several times. Incubate for 5 minutes at 4°C. Centrifuge for 10 seconds.
15. Repeat step 14 again. The two eluates may be pooled in one sampling tube.


(Note: Steps 14-17 of the first protocol are identical to steps 12-15 of the second protocol.)

Related Products:

3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3305A	c-Myc-tagged Protein Mild Purification Kit ver.2	2 purifications
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3307	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 5	
3300-205	c-Myc-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3310A	His-tagged Protein Purification Kit	2 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 2 mg x 5	
3312	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 2 mg x 25	
3310-205	His-tag peptide	2 mg x 5
3317	V5-tagged Protein Purification Kit Ver.2	20 purifications
3317A	V5-tagged Protein Purification Kit Ver.2	2 purifications
3318	V5-tagged Protein Purification Gel Ver.2	1 mL
3315-205	V5-tag peptide	2 mg x 5
3320	HA-tagged Protein Purification Kit	20 purifications
3320A	HA-tagged Protein Purification Kit	2 purifications
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5
3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit	2 purifications
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 25	
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3325-205	DDDDK-tag peptide	1 mg x 5

Example of Purification Results

Purification of c-Myc tagged SIRT1/Sir2 Deacetylase

SDS-PAGE (Coomassie Brilliant Blue Staining)


Human embryonic kidney cells (293T) were transfected with pCMV-c-Myc-SIRT1/Sir2 Deacetylase and cultured for 70 hours. Cells were then lysed in the Lysis buffer (500 μ L/100 mm dish) without any protease inhibitor and purified according to the preceding protocol. For comparison, elution was carried out not only with peptide solution but also with acid solution and alkaline solution. Each purification was conducted with the same amount of Anti-c-Myc tag Beads (5 μ L) and the same amount of elution solution (20 μ L x 2 times and then pooled).

Peptide elution solution: neutral pH

Acid elution solution: pH 2.8 (Neutralize the elution immediately)

Alkali elution solution: pH 11.3 (Neutralize the elution immediately)

Protein Activity

Protein activity was measured by SIRT1/Sir2 Deacetylase Fluorometric Assay Kit.

Additional Information

Several reagents were examined whether or not they were suitable for use with the c-Myc tagged Protein MILD PURIFICATION KIT. The results are listed below.

Chaotropic agents

	Urea	1 M	Yes
	Guanidine-HCl	1 M	No

Reducing agents

	DTT	10 mM	Yes
	2-Mercaptoethanol	10 mM	Yes

Surfactants

Nonionic	Tween-20	5%	Yes
	TritonX-100	5%	Yes
	NP40	5%	Yes
	Digitonin	1%	Yes
	n-Octyl-β-D-gulcoside	1%	Yes
Zwitterionic	CHAPS	1%	Yes
	CHAPSO	1%	Yes
Anionic	SDS	0.05%	No
	Sodium Deoxycholate	0.1%	Yes
	Sodium Deoxycholate	0.5%	No

Others

	NaCl	1 M	Yes
	Glycerol	10%	Yes
	EDTA	10 mM	Yes

The “Yes” indicates the reagents can be used in the Lysis buffer for this kit up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this kit at the indicated concentration.

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