For Research Use Only. Not for use in diagnostic procedures.

# **DDDDK tagged Protein PURIFICATION KIT**

(MoAb. clone FLA-1)

CODE No. 3325

*PURIFICATION to Maintain Protein Activity  
from eukaryote cell lysate and culture supernatant*

## **Product Description**

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more. Thus, choosing a method of purification is an important aspect in maintaining protein structure and function.

Recombinant tagged protein purification methods and kits are now widely recognized. The DDDDK epitope tag (DYKDDDDK) has been widely used as a multi-purpose tag, and anti-DDDDK tag antibodies are optimally suited for identifying, detecting, purifying, and monitoring the expression levels of recombinant DDDDK tagged proteins.

MBL's DDDDK tagged Protein PURIFICATION KIT is designed for the isolation of DDDDK tagged protein from cell culture supernatants and cell lysate under neutral pH condition. Severe conditions such as acidic or alkaline elution denature protein structure. However, a neutral pH elution can preserve protein activity and native conformation. MBL has developed the Anti-DDDDK tag Beads to purify DDDDK tagged proteins quickly and efficiently. As the Beads can be used at neutral pH, the purified proteins can maintain the activity and conformation. The elution of DDDDK tagged proteins from the Beads is achieved by the addition of the DDDDK tag peptide (DYKDDDDK). As the DDDDK tag peptide competes with DDDDK tagged proteins on the Beads, the purified proteins do not lose the protein activity. Using a Spin Column resulting in high efficiency has optimized the simple procedure of this kit.

## **Kit Components**

Components sufficient for conducting 20 times purifications of DDDDK tagged protein.

- |                                |  |
|--------------------------------|--|
| <b>1. Anti-DDDDK tag Beads</b> | 25% slurry: 100 $\mu$ L beads in 400 $\mu$ L total volume in PBS with 0.1% ProClin 150 as preservative |
| <b>2. Elution Peptide</b>      | DDDDK peptide, 1 mg in 1 mL PBS after reconstitution   |
| <b>3. Spin Columns Sets</b>    | 20 columns with pre-inserted bottom plugs and top caps   |
| <b>4. Wash Concentrate</b>     | 10 x concentrate, 6 mL x 2 tubes   |

## **Storage**

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

## **Product Capacity**

The purification capacity of the Anti-DDDDK tag Beads varies depending upon the DDDDK tagged protein.

For examples, 5  $\mu$ L of Anti-DDDDK tag Beads (20  $\mu$ L slurry) bound 10  $\mu$ g of a DDDDK tagged protein (32 kDa) and eluted 9  $\mu$ g of purified protein.

## **Materials Required but not Provided**

1. Micro centrifuge capable of 15,000 x g
2. Sampling tube (1.5 mL)
3. End-over-end rotator
4. Distilled water
5. Lysis buffer  
Suitable Lysis buffer varies with cell type.  
**Note:** see Additional Information

### Homemade Lysis buffer

20-50 mM	Tris-HCl (pH7.5) or HEPES-KOH (pH 7.5)
100-300 mM	NaCl
1-5 mM	EDTA
1%	NP-40 or Triton X-100
	if necessary add Protease Inhibitor Cocktail (e.g. SIGMA: code P8340, PIERCE: code 78415)

## **Protocols**

The following protocols are for the isolation of DDDDK tagged proteins produced in a 100-mm cell culture dish. The expression level of the DDDDK tagged protein may vary. If necessary, adjust the volume of Anti-DDDDK tag Beads and Elution Peptide Solution proportionally.

## **Material Preparation**


### 1. Wash Solution

Dilute Wash Concentrate with 9 times its volume of distilled water.  
(e.g. Dilute 0.5 mL of Wash Concentrate with 4.5 mL of distilled water.)  
For each Spin Column, prepare 5 mL of Wash Solution.

### 2. Elution Peptide Solution

Reconstitute the Elution Peptide with 1 mL of distilled water. If you want to store the reconstituted Elution Peptide, prepare appropriate aliquots (e.g. 45  $\mu$ L x 20 tubes) and store at -20°C. Repeated freezing and thawing is not recommended.

## **Procedure Summary (Purification from mammalian cultured cell lysate)**

## **Purification from mammalian cultured cell lysate**

### **(Lysis of Mammalian Cells)**

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 0.5 mL of Lysis buffer to the cell pellet and vortex.
6. Sonicate the sample for 15 seconds.
7. Incubate the sample for 15 minutes on ice.
8. Remove cell debris by centrifugation at 15,000 x g for 5 minutes at 4°C.

### **(Purification of DDDDK tagged Protein)**


9. Transfer the 0.5 mL of cell lysate (supernatant from step 8) to a 1.5 mL sampling tube.
10. Resuspend the Anti-DDDDK tag Beads by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
11. Dispense 20  $\mu$ L Anti-DDDDK tag Beads suspension (5  $\mu$ L Beads) into the sampling tube (step 9).
12. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
13. Centrifuge at 2,500 x g for 10 seconds. Carefully remove and discard the supernatant (or save for future analysis) using a micropipettor. Don't aspirate beads or disturb bead pellet.
14. Add 1 mL of Wash Solution to the sampling tube. Inverting the sampling tube several times and centrifuge at 2,500 x g for 10 seconds. Carefully remove and discard the supernatant using a micropipettor. Don't aspirate beads or disturb bead pellet. Repeat this step two additional times.
15. Add 0.2 mL of Wash Solution to the sampling tube. Resuspend the Beads by pipetting up and down several times. Transfer the resuspended the Beads to the Spin Column. Take off the bottom plug.
16. Centrifuge at 2,500 x g for 10 seconds. Discard the flow-through.
17. Place the Spin Column in a new sampling tube.
18. For the first elution, screw the bottom plug on tightly. Add 20  $\mu$ L Elution Peptide Solution to the Anti-DDDDK tag Beads, then screw the top cap on tightly. Tap the tube gently several times. Incubate for 5 minutes at 4°C. Remove the top cap and bottom plug. Centrifuge at 2,500 x g for 10 seconds.
19. For the second elution, add 20  $\mu$ L Elution Peptide Solution to the Anti-DDDDK tag Beads, then tap the tube gently several times. Incubate for 1 minute at 4°C. Centrifuge at 2,500 x g for 10 seconds. The two eluates (step 18 and 19) may be pooled in one sampling tube.

### **Related Products:**

3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit	2 purifications
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 25	
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3325-205	DDDDK-tag peptide	1 mg x 5
3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3307	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 5	
3300-205	c-Myc-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 2 mg x 5	
3312	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 2 mg x 25	
3310-205	His-tag peptide	2 mg x 5
3317	V5-tagged Protein Purification Kit Ver.2	20 purifications
3318	V5-tagged Protein Purification Gel Ver.2	1 mL
3315-205	V5-tag peptide	2 mg x 5
3320	HA-tagged Protein Purification Kit	20 purifications
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5

## Example of Purification Results 1

### Purification and enzymatic activity of N-terminal DDDDK tagged $\beta$ -galactosidase

### SDS-PAGE (Coomassie Brilliant Blue Staining) & X-gal staining

Human embryonic kidney cells (293T) were transfected with pcDNA-DDDDK tagged  $\beta$ -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (0.5 mL/100-mm dish) and purified according to the preceding protocol. For comparison, elution was carried out not only with peptide solution but also with acid solution and alkaline solution. Each purification was conducted with the same amount of Anti-DDDDK tag Beads (5  $\mu$ L) and the same amount of elution solution (20  $\mu$ L x 2 times and then pooled).

Peptide elution solution: neutral pH


0.1 M Glycine-HCl: pH 3.0 (Neutralize the eluate immediately with 1 M Tris-HCl, pH 8.0)  
(Acid elution solution)

0.1 M  $\text{NH}_3$ : pH 11.3 (Neutralize the eluate immediately with 1 M acetic acid)  
(Alkali elution solution)

Enzymatic activity of each purification was performed using standard X-gal staining method.

## Example of Purification Results 2

### Purification of Internal DDDDK tagged GFP from *E.coli* lysate

SDS-PAGE (Coomassie Brilliant Blue Staining)

DDDDK tagged GFP plasmid was transformed into *E.coli* BL21(DE3)RIL and cultured in 1 mL LB media at 27°C for 12 hours.

Cells were then lysed in the 0.5 mL Lysis buffer and purified according to the preceding protocol.

## Additional Information

Several reagents were examined whether or not they were suitable for use with the DDDDK tagged Protein PURIFICATION KIT. For example, RIPA buffer could be used for preparation of cell lysate. The results are listed below.

<b>Chaotropic agents</b>			
	<b>Urea</b>	<b>1 M</b>	<b>Yes</b>
	<b>Guanidine-HCl</b>	<b>1 M</b>	<b>No</b>
<b>Reducing agents</b>			
	<b>DTT</b>	<b>10 mM</b>	<b>Yes</b>
	<b>2-Mercaptoethanol</b>	<b>10 mM</b>	<b>Yes</b>
<b>Surfactants</b>			
<b>Nonionic</b>	<b>Tween-20</b>	<b>5%</b>	<b>Yes</b>
	<b>Triton X-100</b>	<b>5%</b>	<b>Yes</b>
	<b>NP-40</b>	<b>1%</b>	<b>Yes</b>
	<b>Digitonin</b>	<b>1%</b>	<b>Yes</b>
	<b>n-Octyl-<math>\beta</math>-D-glucoside</b>	<b>1%</b>	<b>Yes</b>
<b>Zwitterionic</b>	<b>CHAPS</b>	<b>1%</b>	<b>Yes</b>
	<b>CHAPSO</b>	<b>1%</b>	<b>Yes</b>
<b>Anionic</b>	<b>SDS</b>	<b>0.1%</b>	<b>Yes</b>
	<b>Sodium Deoxycholate</b>	<b>0.5%</b>	<b>Yes</b>
<b>Others</b>			
	<b>NaCl</b>	<b>1 M</b>	<b>Yes</b>
	<b>Glycerol</b>	<b>10%</b>	<b>Yes</b>
	<b>EDTA</b>	<b>10 mM</b>	<b>Yes</b>

The “Yes” indicates the reagents can be used in the Lysis buffer for this kit up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this kit at the indicated concentration.

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3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit	2 purifications
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 1, Peptide: 1 mg x 5
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 5, Peptide: 1 mg x 25
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3325-205	DDDDK-tag peptide	1 mg x 5
3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	Gel: 1 mL x 1, Peptide: 1 mg x 5
3307	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	Gel: 1 mL x 5, Peptide: 1 mg x 5
3300-205	c-Myc-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 1, Peptide: 2 mg x 5
3312	His-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 5, Peptide: 2 mg x 25
3310-205	His-tag peptide	2 mg x 5
3317	V5-tagged Protein Purification Kit Ver.2	20 purifications
3318	V5-tagged Protein Purification Gel Ver.2	1 mL
3315-205	V5-tag peptide	2 mg x 5
3320	HA-tagged Protein Purification Kit	20 purifications
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5

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