CODE No. 4820

MBL

Printed: November 15, 2004

Revised: September 1, 2011

For Research Use Only. Not for use in diagnostic procedures.

Detection of caspase activity

APOPCYTO

Caspase-8 Fluorometric Assay Kit

CODE No. 4820

100 assays

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
URL <https://ruo.mbl.co.jp>
e-mail support@mbi.co.jp, TEL 052-238-1904

CODE No. 4820

Before use, thoroughly read these Instructions.

Description

Caspase is a member of the cysteine aspartic acid-specific protease family, which is activated by a variety of signals of death receptor ligation, DNA damages, serum starvation and stresses *etc.* Active caspase recognizes a lot of several molecules as substrates to cleave them, occurring to biological events corresponding to the apoptosis. For example, ICAD (inhibitor of caspase-activated deoxyribonuclease) is inactivated and CAD (caspase-activated deoxyribonuclease) is indirectly activated by caspase-3, and it is related to chromatin fragmentation for nucleosome units. Caspase recognizes several structural proteins as a substrate to cleave them, and the cleavage is associated with the unique apoptosis cell morphology of chromatin condensation, nucleus fragmentation and cytoplasmic integrity. Active caspase could act in apoptosis process, so that detection of caspase activity is necessary for apoptosis research frequently. Though there are numerous ways to measure of caspase activity, colorimetric or fluorometric substrates are used popularly.

Assay Principle

Active caspase recognizes 4 amino acid residues in a substrate molecule, which are from the aspartate residue next to N-terminal 3 amino acid residues, and specifically cleaves at the C-terminal side of the aspartate residue. So this characterization is useful for detecting caspase activity by using synthetic substrates. In this method, 4 amino acid sequences is labeled with AMC (7-amino-4-methyl-coumarin) pNA (*p*-nitroanilide), or AFC (7-amino-4-trifluoromethyl coumarin) at the C-terminal side. Free pNA, AMC or AFC is released from the labeled synthetic substrate on cleavage by active caspase, and monitored by a 96-well microplate reader (pNA: absorbance wavelength 400 or 405 nm, AMC: excitation wavelength 380 nm, fluorescence wavelength 460 nm, AFC: excitation wavelength 400 nm, fluorescence wavelength 505 nm). For example, if caspase is active in sample, tetrapeptide-AMC as a synthetic substrate is cleaved and the excitation wavelength and the fluorescence wavelength are shifted. Because the amount of free AMC is proportional to the amount of caspase activity present in the sample, caspase activity can be calculated by monitoring of the optical density of free AMC at wavelength Excitation 380 nm/Emission 460 nm.

APOPCYTO Caspase-8 Fluorometric Assay Kit is provided to detect caspase-8 activity in cell extract with IETD-AMC as a substrate, which IETD sequence is recognized by active caspase-8 selectively. APOPCYTO Caspase-8 Fluorometric Assay Kit is composed of reagents necessary for the detection, which are Cell Lysis buffer, Reaction buffer and so on. Caspase-8 activity can be detected easily and quickly using APOPCYTO Caspase-8 Fluorometric Assay Kit.

Intended Use

For research use only. Not for use in *in vitro* diagnostic procedures for clinical diagnosis.

Product Components

Materials	Quantity	Color
Cell Lysis Buffer	50 mL × 1 vial	Clear
2X Reaction Buffer	2 mL × 3 vials	Blue
1 M DTT	500 µL × 1 vial	Yellow
Substrate IETD-AMC (2 mM)	550 µL × 1 vial	Brown
Inhibitor IETD-FMK (1 mM)	55 µL × 1 vial	Orange
10 mM AMC	200 µL × 1 vial	Green

Expiration

Please see the label of this kit.

Storage Conditions

Store at -20°C. After thaw, the Cell Lysis Buffer and 2X Reaction Buffer can be stored at 4°C.

Material to Be Supplied by the User

- microcentrifuge (For harvest cells)
- 96-well microplate reader
(For measure the fluorescence counts, Excitation 380 nm/Emission 460 nm)
- micropipette
- 96-well microplate (flat-bottom, clear polystyrene)
- microcentrifuge tube

Note

- When caspase-8 activity is not measured immediately after cell preparations, store the cell lysates or cell pellets at -20°C and use them within 1 week.
- Store caspase-8 substrate away from light.
- Ensure that DTT is added to 2X Reaction Buffer to obtain a final concentration of 10 mM before use.
- We recommend you to previously measure the fluorescence counts in the empty wells of 96-well microplate with a 380 nm excitation filter and 460 nm emission filter to avoid use of any wells with high background.

Preparation of reagents

- 1) Add DTT to 2X Reaction Buffer to obtain final concentration of 10 mM just before use. (For example, 10 μ L of 1 M DTT to 990 μ L of 2X Reaction Buffer)
- 2) Preparation of AMC (7-amino-4-methyl-coumarin) standards
 - i) Dilute 10 mM AMC to 500 μ M AMC in Cell Lysis Buffer (5 μ L of 10 mM AMC is added to 95 μ L of Cell Lysis Buffer).
 - ii) Dilute 500 μ M AMC in Cell Lysis Buffer according to the following table to give these final concentrations.

Final AMC concentration	AMC	Cell Lysis Buffer
50 μ M (5 nmol)	50 μ L of 500 μ M AMC	450 μ L
25 μ M (2.5 nmol)	250 μ L of 50 μ M AMC	250 μ L
12.5 μ M (1.25 nmol)	250 μ L of 25 μ M AMC	250 μ L
6.25 μ M (0.625 nmol)	250 μ L of 12.5 μ M AMC	250 μ L
3.125 μ M (0.312 nmol)	250 μ L of 6.25 μ M AMC	250 μ L
1.5625 μ M (0.165 nmol)	250 μ L of 3.125 μ M AMC	250 μ L
0 μ M (0 nmol)		250 μ L

Caspase Assay Protocol

1. Induce apoptosis in cells using the desired method, Concurrently incubate a control culture without induction.
2. Count the cells and harvest the $1-5 \times 10^6$ cells by centrifugation at 400 X *g* for 5 minutes.
3. Resuspend the cells in ice-cold Cell Lysis Buffer (50-500 μ L), and incubate the cells on ice for 10 minutes.
4. Centrifuge the cell lysates at 10,000 X *g* for 5 minutes at 4°C to precipitate cellular debris. Transfer the supernatants (cell extracts) to new microcentrifuge tubes and put on ice.

If necessary, total protein concentration of the cell extracts may be measured by standard protein assay methods, and adjust protein concentration to 100-200 μ g/50 μ L in Cell Lysis Buffer.

5. Add 50 μ L of 2X Reaction Buffer containing 10 mM DTT to each well of a 96-well microplate.
10 mM DTT should be added to 2X Reaction Buffer just before use.
6. Add 45 μ L of cell lysates to each well. Prepare wells which were added 45 μ L of Cell Lysis Buffer in stead of cell lysates to measure blank absorbance.

[Optional method]

To verify that the signal detected by the kit is due to protease activity, incubate an induced sample with IETD-FMK before adding substrate. Add 45 μ L of cell extracts to 50 μ L of 2X Reaction Buffer (containing 10 mM DTT) and 1 μ L of 1 mM IETD-FMK.

8. Add 5 μL of Caspase-8 Substrate to each well.
9. Cover the plate with a plate sealer and incubate at 37 $^{\circ}\text{C}$ for 60 minutes.
10. Add 100 μL of AMC standards to empty wells.
11. Measure the fluorescence counts in the wells with 380 nm excitation filter and 460 nm emission filter. Calculate the specific activity (SA) of caspase-8 present in each sample using following formula:
 - i) Construct a standard curve using the fluorescence counts of AMC (7-amino-4-methyl-coumarin) standards.
 - ii) Calculate the relative absorbance (Fr) as:

$$\text{Fa} = (\text{induced apoptosis sample counts}) - (\text{blank counts})$$

$$\text{Fn} = (\text{negative sample counts}) - (\text{blank counts})$$

$$\text{Fr} = \text{Fa} - \text{Fn}$$
 - iii) Calculate liberated free AMC concentration (B μM) equal to Fr by AMC standard curve.
 - iv) Calculate caspase-8 activity (C) present in each sample as:

$$\text{C} = \text{nmole AMC liberated per hour}$$

$$= \text{BX } 0.1 \text{ mL} / \text{incubation hour}$$

$$= 0.1\text{B} / \text{incubation hour}$$
 - v) If total protein concentration of the cell extracts is D mg/mL, calculate the specific activity (SA) of caspase-8 present in each sample as:

$$\text{SA} = \text{C per } \mu\text{g protein}$$

$$= (0.1\text{B} / \text{incubation hour}) / (\text{D mg/mL} \times 0.1 \text{ mL})$$

$$= \text{B} / (\text{D} \cdot \text{incubation hour})$$

Assay Example


Fig.1 Caspase-8 activity in CH-11 treated Jurkat cells

After Jurkat cells were incubated with 100 ng/mL of anti-Fas monoclonal antibody (CH-11) at 37 $^{\circ}\text{C}$ for 4 hours. Caspase-8 activity was measured by APOPCYTO Caspase-8 Fluorometric Assay Kit.

References

- 1) Wolf, B. B., *et al.*, *J. Biol. Chem.* **274**, 20049-20052 (1999)
- 2) Thornberry, N. A., *et al.*, *J. Biol. Chem.* **272**, 17907-17911 (1997)
- 3) Walker, N. P., *et al.*, *Cell* **78**, 343-352 (1994)
- 4) Sleath, P. R., *et al.*, *J. Biol. Chem.* **265**, 14526-14528 (1990)

Manufacturer

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL <https://ruo.mbl.co.jp>

e-mail support@mbi.co.jp, TEL 052-238-1904

Before use, thoroughly read these Instructions.

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