Printed: July, 2016

For Research Use Only. Not for use in diagnostic procedures.

In situ detection kit for Programmed Cell Death

MEBSTAIN
Apoptosis TUNEL Kit III

CODE No. 8442

FCM: 66 tests/stain: 40 tests

CONTENTS

1. Intended Use -----	1
2. Summary and Explanation-----	1
3. Principle -----	1
4. Materials provided -----	2
5. Storage and Stability -----	2
6. Precautions -----	2
7. Procedure -----	2
8. Related products-----	6
9. References -----	7

Before use, thoroughly read these Instructions.

Intended Use

The MEBSTAIN[®] Apoptosis Kit III is intended for the detection of apoptotic cells in research samples. It is suitable for histochemistry, cytospin preparation, and flow cytometry. **For Research Use Only. Not for use in diagnostic procedures.**


Summary and Explanation

Programmed cell death is a selective process of physiological cell deletion. It plays a major role in developmental biology and in maintenance of homeostasis in vertebrates^{8), 9), 12), 13)}. For example, during maturation, T cells recognizing self antigens are destroyed by programmed cell death¹¹⁾. Programmed cell death is morphologically known as apoptosis. Apoptosis is accompanied by condensation of cytoplasm, loss of plasma membrane microvilli, condensation and fragmentation of nuclei, and extensive degradation of chromosomal DNA into oligomers of about 180 bp. It has been considered that fragmentation of nuclear DNA is a biochemical hallmark of apoptosis. The TdT-mediated dUTP-biotin nick end labeling method (TUNEL method) developed by Gavrieli et al.¹⁴⁾ has enabled *in situ* visualization of DNA fragmentation at the single cell level and is considered to be a more sensitive method than conventional morphological techniques¹⁰⁾.

The MEBSTAIN Apoptosis Kit III is an apoptosis detection kit based on the TUNEL method.

Principle

In cells in which apoptosis occurs, the chromatin DNA is cut by endonuclease at linker DNA site between nucleosomes. Then, DNA fragments, which are a number of multimers of nucleosomal units exit in nuclei. In TUNEL method, 3'-OH DNA ends generated by DNA fragmentation is nick end labeled with Biotin-dUTP, mediated by terminal deoxynucleotidyl transferase (TdT), and then stained using Avidin-DTAF to allow specific staining.

Materials provided

Each kit contains;

Name	Materials	Quantity (for 1 kit)
Proteinase K (PK)	Proteinase K	0.4 mL × 1 vial
TdT	TdT	0.1 mL × 1 vial
TdT buffer II	Buffer for preparing TdT solution	2.0 mL × 3 vials
Biotin-dUTP II	Biotin conjugated dUTP (ready-for-use)	0.1 mL × 1 vial
TB buffer (×10)	TB buffer (×10)	50 mL × 1 vial
Blocking concentrate (×4)	Buffer for blocking samples (×4)	2.0 mL × 1 vial
Avidin-DTAF	DTAF conjugated streptavidin	0.2 mL × 1 vial

Storage and Stability

All kit components must be stored at -20°C. All reagents are stable until the expiration date when stored at the indicated conditions.

Precautions

1. Allow all the components to come to room temperature (20-25°C) before use.
2. **Blocking concentrate** and **Avidin-DTAF** contain sodium azide (<0.1%) as preservative. Azide may react with copper or lead in plumbing system to form explosive metal azide. Therefore, always flush with plenty of water into a drain when disposing materials containing azide.
3. **TdT** contains 90 mM sodium cacodylate and 90 mM potassium cacodylate. **TdT buffer II** contains 100 mM potassium cacodylate. As cacodylic acid is arsenic compound, do not dispose materials containing TdT and TdT buffer II into a drain.
4. **TdT** contains 2-mercaptoethanol (0.05%). Handle with care.
5. This kit is intended for research use only. Not for use in diagnostic procedures.

Procedure

◆ Preparation of Reagents

1. Proteinase K (PK) solution

Prepare “**PK solution**” by diluting **Proteinase K, 1:25** with PBS. Warm this solution to 37°C prior to use.

Example: Dilute 10 µL of Proteinase K to 240 µL of PBS.

2. TdT solution

Prepare “**TdT solution**” by mixing the **TdT buffer II, Biotin-dUTP II, and TdT** at **18:1:1** prior to use.

CODE No. 8442

Example: Mix 45 μ L of TdT buffer II, 2.5 μ L of Biotin-dUTP II, and 2.5 μ L of TdT.

* Prepare “**TdT solution**” on the ice.

3. TB solution

Prepare “**TB solution**” by diluting **TB buffer, 1:10** with distilled water prior to use.

Example: Dilute 50 mL of TB buffer to 450 mL of distilled water.

* After thawing TB buffer, it should be stored at 4°C.

4. Blocking solution

Prepare “**Blocking solution**” by diluting **Blocking concentrate, 1:4** with PBS prior to use.

Example: Dilute 25 μ L of Blocking concentrate to 75 μ L of PBS.

* After thawing Blocking concentrate, it should be stored at 4°C.

5. Avidin-DTAF solution

Prepare “**Avidin-DTAF solution**” by diluting **Avidin-DTAF, 1:10** with Blocking solution prior to use.

Example: Dilute 10 μ L of Avidin-DTAF to 90 μ L of Blocking solution.

* After thawing Blocking concentrate, it should be stored at 4°C.

◆ for Immunohistochemistry

Materials required but not provided.

- Xylen
- Ethanol (100%, 90%, 80%)
- Phosphate buffered saline (PBS)
- Distilled water
- Mounting medium (90% glycerin, 10% PBS)
- 4% paraformaldehyde [0.1M NaH_2PO_4 , pH7.4]
- Propidium Iodide (0.5 mg/mL) *It is for counterstaining.
- Cover slip
- Moist chamber
- Adequate test tube
- Adjustable micropipette

STEP 1. (Preparation of tissue sections)

- 1) Use glass slides pretreated with silane.
- 2) Fix the samples with 4% paraformaldehyde.

Note

4% paraformaldehyde is recommended for sensitive detection. Do not use ethanol and acetone because they may cause non-specific staining.

3-5 μm thick of sections are suitable for staining.

a) Paraffin embedded tissue sections

Fix tissue samples in 4% paraformaldehyde and embedded in paraffin.

b) Frozen tissue sections

Fix cryosections in 4% paraformaldehyde.

When using frozen tissue sections, be careful not to allow it to dry. Proceed directly **STEP 4.**

(DNA nick end labeling).

STEP 2. (Deparaffinization)

Be careful not to allow tissue sections to dry at each step.

- 1) Deparaffinize the sections with Xylene 3 times for 3 minutes each.
- 2) Wash the slides with 100% Ethanol 3 times for 3 minutes each.
- 3) Wash the slides with 90% Ethanol 3 times for 3 minutes each.
- 4) Wash the slides with 80% Ethanol 3 times for 3 minutes each.
- 5) Wash the slides with distilled water 4 times for 2 minutes each.

STEP 3. (Treatment with proteinase K (PK))

This step is not required for frozen tissue sections.

- 1) Pipette 500 μL of PBS on the section and incubate it for 30 minutes at 37°C. (Or put the slide into a jar filled with PBS which is warmed to 37°C, and incubate it for 15-30 minutes at 37°C.)
- 2) After removing PBS buffer by tapping off, pipette 250 μL of **PK solution** on the section and incubate it for 30 minutes at 37°C.
- 3) Wash the slides with distilled water at 4 times for 2 minutes each.

Note

In the meantime, prepare **TdT solution**. Preparation should be in the ice.

STEP 4. (DNA nick end labeling)

- 1) Pipette 50 μL of **TdT buffer II** on the section and incubate it for 5-10 minutes at room temperature.
- 2) After removing **TdT buffer II** by tapping off, pipette 50 μL of **TdT solution** and incubate it for 60 minutes at 37°C.

CODE No. 8442

- 3) After removing **TdT solution** by tapping off, immerse the slides in **TB solution** in a coplin jar and incubate it for 15 minutes at room temperature.
- 4) Wash the slides with distilled water 4 times for 2 minutes each.

STEP 5. (Histchemistry)

- 1) Pipette 50 μL of **Blocking solution** on the section and incubate it for 10 minutes at room temperature.
- 2) After removing **Blocking solution** by tapping off, pipette 50 μL of **Avidin-DTAF solution** on the section and incubate it for 30 minutes at 37°C.
- 3) Wash the slides with PBS at 3 times for 5 minutes each at room temperature.

STEP 6. (Counterstaining)

- 1) Pipette 50 μL of Propidium Iodide (0.5 $\mu\text{g}/\text{mL}$) on the section and incubate it for 15-20 minutes at 4°C.
- 2) Wash the slides with PBS for 5 minutes.

STEP 7. (Mounting and microscopy)

- 1) Mount with mounting medium (90% glycerin, 10% PBS) under coverslip.
- 2) View by fluorescent microscopy.

◆ for Flow Cytometry

It is recommended that the positive and negative controls are set.

Example: as a negative control; peripheral blood mono nuclear cell or cultured cells, such as Raji cells etc. (more than 95% viability).

as a positive control; those above cells which are treated with DNase I (1 $\mu\text{g}/\text{mL}$) for 60 minutes at 37°C. DNase I treatment should be done after permeabilization.

Materials required but not provided.

- Phosphate buffered saline (PBS)
- PBS containing 0.2% BSA
- 4% paraformaldehyde [0.1M NaH_2PO_4 , pH7.4]
- 70% Ethanol
- Distilled water
- Adequate test tube
- Adjustable micropipette

STEP 1. (Fixation of cells)

- 1) Wash the cells with PBS containing 0.2% BSA several times.
- 2) Fix the cells with 4% paraformaldehyde for 30 minutes at 4°C.
- 3) Wash the cells with PBS containing 0.2% BSA 2 times.

Note

4% paraformaldehyde is recommended for sensitive detection. Do not use ethanol and acetone because they may cause non-specific staining.

STEP 2. (Permeabilization)

Add 200 µL of 70% Ethanol to the cell pellet, and mix it gently. Then incubate it for 30 minutes at -20°C to permeabilize.

Note

In the meantime, prepare TdT solution. Preparation should be done in the ice.

When the next step can't be done right after this step, cells are suspended in 70% Ethanol and stored at -20°C until test. The test should be done within one week.

STEP 3. (DNA nick end labeling)

- 1) After washing with PBS containing 0.2% BSA 2 times, add 30 µL of **TdT solution** to the cell pellet and incubate it for 60 minutes at 37°C. As negative control, add 30 µL of **negative TdT solution (TdT buffer II and Biotin-dUTP II at 19:1)** to the cell pellet and incubate it for 60 minutes at 37°C.
- 2) Wash the cells with PBS containing 0.2% BSA 2 times.

STEP 4. (Staining)

- 1) Add 10 µL of **Blocking solution** to the cell pellet then keep it for 10 minutes at room temperature.
- 2) Add 30 µL of **Avidin-DTAF solution**, then incubate it for 30 minutes at room temperature.
- 3) Wash the cells with PBS containing 0.2% BSA 2 times.
- 4) Suspend the cell pellet to 200-500 µL of PBS containing 0.2% BSA and analyze by a flow cytometer.

Related products

8445	MEBSTAIN Apoptosis TUNEL Kit Direct
4700	MEBCYTO [®] Apoptosis Kit (Annexin V-FITC Kit)
SY-001	Anti-Fas (CD95) mAb (clone CH-11)

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