Anti-His-tag mAb-HRP-Direct

Code No.	Clone	Subclass	Quantity
D291-7	OGHis	Mouse IgG2aκ	50 μL

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HRP-Direct Series

HRP-Direct is a series of HRP conjugated primary antibodies developed by MBL.

HRP-Direct Series products don't need secondary antibodies. That brings the following advantages:

- ① Total incubation time is cut in half.

The reaction with the secondary antibody becomes unnecessary. In the shaking method, reaction time is only 30 minutes. Spend your evening doing more important things.

- ② Clear result

No more heavy and light chain bands in Immunoprecipitation. The HRP conjugated primary antibody will not detect your precipitating antibody. Clear result helps your research.

Usually, the secondary antibody amplifies the signal of an unconjugated primary antibody. Therefore, it is thought to be difficult to obtain a strong signal with a directly conjugated antibody. To overcome this hurdle, MBL has improved the HRP conjugation method. MBL's HRP-Direct Series yield strong signal with minimum background. Give it a try - the HRP-Direct Series will not disappoint you.

DETECTED ANTIGEN

N-terminal, Internal, C-terminal His-tag

REAGENT

PBS/Preservative/Stabilizer

STORAGE

This antibody is stable for one year from the date of purchase when stored at 4°C.

INTENDED USE

For Research Use Only. Not for use in diagnostic procedures.


APPLICATION

Western blotting;

Method	Reaction times	Dilution	Detection
Sitting	1 hour	5,000	chemiluminescence
Shaking	30 minutes	5,000	chemiluminescence


*Detailed procedure is provided in the following protocols.

- Comparison of Sitting method and Shaking method

Lane1: N-terminal Met-His-tagged protein A expressed in *E. coli*
 Lane2: Internal His-tagged protein B expressed in 293T
 Lane3: Internal His-tagged protein C expressed in *E. coli*
 Lane4: C-terminal His-tagged protein D expressed in 293T
 Lane5: C-terminal His-tagged protein E expressed in 293T

- IP-Western

Immunocytochemistry; Not tested

Immunohistochemistry; Not tested

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

PROTOCOL

Western blotting

1. SDS-PAGE and Western blotting

Perform SDS-PAGE electrophoresis on the protein samples and transfer the protein to a PVDF membrane according to standard techniques.

2. Blocking

Block the membrane with 5% skimmed milk in PBS (or TBS). Incubate for 1 hour with shaking at room temperature or overnight at 4°C.


3. Incubate the membrane with HRP-Direct antibody as follows:

① Sitting Method (Saving type: amount of antibody)

Dilute the HRP-Direct antibody with 1% skimmed milk in PBS (or TBS). Use at least 30 $\mu\text{L}/\text{cm}^2$ membrane. Antibody dilution ratio is suggested in the APPLICATION. Drain the excess blocking buffer from the membranes and place them, protein side up, on a wrap film or other suitable clean surface. Pipette the HRP-Direct solution on to the membrane. Cover the membrane with plastic case to avoid drying. Incubate for 1 hour at room temperature.

② Shaking Method (Saving type: reaction time)

Dilute the HRP-Direct antibody with 1% skimmed milk in PBS (or TBS). Use at least 10 mL. Antibody dilution ratio is suggested in the APPLICATION. Transfer the membrane to a tray containing HRP-Direct antibody. Incubate for 30 minutes with shaking at room temperature.

4. Washing

Wash the membrane in PBS (or TBS) with 0.05% Tween-20 three times for 5 minutes each.

5. Detection

Detect with a chemiluminescence substrate according to the manufacturer's instructions.