MJS002V2 Protein G-Magnetic Beads

PRODUCT DESCRIPTION

Protein G-Magnetic Beads are well-designed magnetic microparticles for bioseparation such as immunoprecipitation and antibody purification. Their surfaces are covered with a JSR Life Sciences proprietary hydrophilic polymer and chemically conjugated Protein G, which specifically binds to the F_C moiety of Immunoglobulin. These chemistries enable you to perform a high yield and low non-specific binding bioseparation in a simple way.

FEATURES

- Uniform particle size
- Superparamagnetic
- Rapid magnetic responsiveness
- High IgG capacity
- Low non-specific binding

EXAMPLE APPLICATIONS

Immunoprecipitation(Protein, Chromatin), Antibody purification

SPECIFICATIONS

Package volume	10 mL
Solid content in slurry	1% (6 x 10 ⁸ beads/mL approx.)
Dispersion media	TBS* + 0.05% Tween20 + 0.09 % Sodium Azide
Bead diameter	3 μm (micrometer)
Bead magnetite content	20% approx.
IgG capacity	7 μg/mg bead approx. (In case of Mouse IgG)
Shelf life	Labeled on the bottle
*TBS: Tris buffered saline,	25 mM Tris-HCl (pH 7.2)/0.15 M NaCl

STORAGE

Protein G-Magnetic Beads are stable for 24 months when stored at 2-8 °C. Do not freeze the vial. Vortex the vial or pipette gently up and down to obtain a homogeneous dispersion before use

DISPOSAL

Reagent contains sodium azide at a low concentration as a preservative. Sodium azide is toxic if ingested and may react with heavy metals to form explosive metal azides. Azide compounds should be diluted with running water before discarding to avoid deposits in plumbing where explosive condition may develop.

RECOMMENDED PROTOCOLS

Protocol I & II for antibody binding onto **Protein G-Magnetic Beads**. (Preparation for Immunoprecipitation)

Protocol-1 (w/o IgG Cross-linking) may have an advantage in terms of maintaining antibody activity.

In case you need to avoid co-elution of your antibody, please try Protocol-2 (Cross-linking). The optimal condition may depend on the nature of your antibody.

[Protocol I] IgG Binding without Cross-linking

Reagent and equipment requirement

Binding Buffer	: Citrate-Phosphate buffer (pH5.0) with 0.1% Tween 20
Washing & Storage Buffer	: TBS-T (25 mM Tris-HCl, pH7.2, 0.15 M NaCl, 0.05% Tween 20)
Equipment	: Magnetic separator. Vortex tube mixer. Tube rotator.

1. Suspend the **Protein G-Magnetic Beads** well using Vortex mixer and put 1 mL of the suspension (i.e., 10 mg beads) into a microtube.
2. Place the tube on a magnetic separator for 1 minute (or longer if needed) and remove the supernatant carefully.
3. Add 1 mL of Binding Buffer and suspend the beads by vortexing. Then, remove the supernatant as in step 2.
4. Add 100 μL of Binding Buffer and 100 μg of IgG solution (100 μL, if antibody was diluted to 1 mg/mL) and suspend the beads by vortexing.
5. Keep rotating the tube using a tube rotator for 30 minutes at room temperature.
6. Remove the supernatant as in step 2.
7. Add 1 mL of Washing Buffer and suspend them by vortexing.
8. Repeat steps 6 & 7 again. Then, remove the supernatant as in step 2.
9. Suspend the beads with Washing & Storage Buffer or a desired buffer suitable for downstream applications and store at 2-8°C until needed.

[Protocol II] IgG Binding with Cross-linking

Reagent and equipment requirement

Binding Buffer	: Citrate-Phosphate buffer (pH 5.0) with 0.1 % Tween 20
Cross-link Buffer	: 0.2 M Triethanolamine (pH 8.2, Adjusted with HCl(aq.))
Cross-link Reagent	: 20 mM DMP; prepared just before the cross-linking reaction (DMP: Dimethyl pimelimidate-2HCl, Thermo Fisher Scientific Inc.)
Stop Solution	: TBS (25 mM Tris + 0.15 M NaCl, pH 7.2)
Washing & Storage Buffer	: TBS-T (25 mM Tris-HCl, pH 7.2, 0.15 M NaCl, 0.05% Tween 20)
Equipment	: Magnetic separator. Vortex tube mixer. Tube rotator.

1. Follow the same procedure found in Protocol I steps 1~6.
2. Add 1 mL of Cross-link Buffer and suspend the beads by vortexing. Then, place the tube on a magnetic separator for 1 minute (or longer if needed) and remove the supernatant carefully. Repeat this step twice.
3. Add 1 mL of Cross-link Reagent and suspend the beads by vortexing.
4. Keep rotating the tube using a tube rotator for 30 minutes at room temperature.
5. Place the tube on magnetic separator for 1 minute and remove the supernatant carefully.
6. Add 1 mL of Stop Solution and suspend the beads by vortexing.
7. Keep rotating the tube using a tube rotator for 30 minutes at room temperature.

8. Remove the supernatant as in step 5.
9. Add 1 mL of Washing Buffer and suspend them by vortexing.
10. Repeat steps 8 & 9 again. Then remove the supernatant as in step 5.
11. Suspend the beads with Washing & Storage Buffer or a desired buffer suitable for downstream applications and store at 2-8 °C until needed.

Protocol III for Immunoprecipitation (IP)

Antibody conjugated beads prepared from Protocol I and Protocol II above are applicable to the following immunoprecipitation protocol.

[Protocol III] Immunoprecipitation of Target Protein

Reagent and equipment requirement

Specimen	: Cell lysate, Serum (10-5,000 µg-protein/mL)
IP Washing Buffer	: 20 mM HEPES(pH 7.9) + 10 v/v% glycerol + 0.5 M KCl + 0.1 % NP-40 + 0.1 mM EDTA
Eluent	: 0.5% Sodium dodecyl sulfate or a desired acid solution suitable for downstream applications
Equipment	: Magnetic separator. Vortex tube mixer. Microtube shaker.

1. Suspend the antibody conjugated **Protein G-Magnetic Beads** well using Vortex mixer and put the 1mg beads (100 µL, if the beads were stored at 10 mg/mL) into a microtube.
2. Place the tube on magnetic separator for 1 minute and remove the supernatant carefully.
3. Add your specimen (20~1000 µL) and suspend the beads by vortexing.
4. Incubate the beads for 30 minutes at room temperature with gentle mixing.
5. Remove the supernatant as in step 2.
6. Add 1 mL of IP Washing Buffer and suspend the beads by vortexing.
7. Repeat steps 5 & 6 for a total of 3 times.
8. Transfer the suspension to a clean tube.
9. Remove the supernatant as in step 2.
10. Add 20 µL of Eluent and suspend the beads by vortexing.
11. Incubate the beads for 10 minutes at room temperature with gentle mixing.
12. Place the tube on a magnetic separator for 1 minute and transfer the eluate to a clean tube.

[Protocol IV] antibody purification

Antibodies like Human IgG, Mouse IgG, Rabbit IgG and so on, which specifically bind to Protein G, can be extracted from your specimen by the following protocol.

Reagent and equipment requirement

Specimen	: Serum
Washing Buffer-1	: Citrate buffer (pH 5.0) + 0.01% Tween 20
Eluent	: Citric acid (pH 2.0)
Washing Buffer-2	: 0.5 M Phosphate buffer (pH 8.0)
Equipment	: Magnetic separator. Vortex tube mixer. Microtube shaker.

1. Suspend the **Protein G-Magnetic Beads** well using a Vortex mixer and put 100 µL of the suspension (i.e., 1 mg beads) into a microtube.
2. Place the tube on magnetic separator for 1 minute and remove the supernatant carefully.
3. Add 1 mL of Washing Buffer-1 and suspend the beads by vortexing.
4. Remove the supernatant as in step 2.
5. Repeat steps 3 & 4 again.
6. Add 200 µL of Specimen and suspend the beads by vortexing.
7. Keep shaking the tube using a Microtube shaker for 20 minutes at room temperature.
8. Remove the supernatant as in step 2.
9. Add 1 mL of Washing Buffer-1 and suspend them by vortexing.
10. Repeat steps 8 & 9 again. Then remove the supernatant as in step 2.
11. Add 100 µL of Eluent and suspend the beads by vortexing.
12. Keep shaking the tube using a Microtube shaker for 20 minutes at room temperature.
13. Place the tube on a magnetic separator for 1 minute and transfer the eluted protein to a clean tube.

In some cases, Protein G beads can be regenerated and reused by the following steps;


14. Wash the beads using 1 mL of Washing Buffer-1 and remove the supernatant.
15. Wash the beads using 1 mL of Washing Buffer-2 and remove the supernatant.
16. Suspend the beads with 1 mL of Washing Buffer-1 and return to the step 3.

EXPERIMENTAL EXAMPLE

[Example I] Immunoprecipitation

Immunoprecipitation of 20S proteasome complex from Jurkat cell lysate
Anti-20S proteasome alfa6 monoclonal antibody (Biomol International, L.P, Clone MCP20) was coupled onto **Protein G-Magnetic Beads** through **Protocol II**.

Immunoprecipitation was performed through **Protocol III**.

Lane 1 Molecular weight marker
Lane 2 IP product using **Protein G-Magnetic Beads**
Lane 3 IP product using competitor-A's magnetic beads
Lane 4 IP product using competitor-B's magnetic beads


<Conditions>

Specimen: Jurkat cell lysate 100µL (30µg protein)
Elution: 20 µL of 0.5% Sodium dodecyl sulfate
Detection: SDS-PAGE (reduced), Silver stain

By using **Protein G-Magnetic Beads**, subunits (alfa1- alfa7, beta1- beta7) of 20S proteasome complex were isolated with high purity from the cell lysate (lane 2, *inside dotted frame*). In comparison, competitor's Protein G magnetic beads resulted in low yield (lane 3) or low purity (lane 4) under the same experimental condition.

[Example II] Antibody purification

Human IgG was extracted from human serum by using **Protein G-Magnetic Beads** through **Protocol IV**:

By using **Protein G-Magnetic Beads**, human antibody was isolated with high purity from the human serum (Lane 2). Recycled beads also showed almost the same performance (Lane 3, 4).

IMPORTANT NOTICE

This product is for research use only and not intended for therapeutic or *in vivo* diagnostic use.
The specificity of this product may change without notice.
JSR Life Sciences Corporation does not guarantee that this product will be continuously available.
JSR Life Sciences Corporation makes no warranties as to this product including, but not limited to, implied warranties of merchantability or fitness for a particular purpose.

This product has used Protein G-Magnetic Beads manufactured by JSR Corporation. Magnosphere™ is a trademark of JSR Corporation.

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