

For Research Use Only.
Not for use in diagnostic procedures.



Anti-Atg2A pAb

CODE No. PD041

CLONALITY Polyclonal
ISOTYPE Rabbit Ig, affinity purified
QUANTITY 100 µL

SOURCE Purified Ig from rabbit serum
IMMUNOGEN Human Atg2A, 700-1,400 aa (recombinant)
FORMULATION PBS containing 50% Glycerol (pH 7.2). No preservative is contained.
STORAGE This antibody solution is stable for one year from the date of purchase when stored at -20°C.

APPLICATIONS-CONFIRMED

Western blotting 1:1,000 for chemiluminescence detection system
Immunoprecipitation 5 µL/300 µL of cell extract from 3 x 10⁶ cells
Immunocytochemistry 1:400

SPECIES CROSS REACTIVITY on WB

Species	Human	Mouse	Rat	Hamster
Cells	293T	MEF	PC12	CHO
Reactivity	+	+	+	+

Entrez Gene ID 23130 (Human), 329015 (Mouse), 689688 (Rat), 100758391 (Hamster)

REFERENCES

- 1) Tamura, N., *et al.* *FEBS Lett.* **591**, 3819-3830 (2017) [WB]
- 2) Tang, Z., *et al.*, *Cell Death Differ.* **24**, 2127-2138 (2017) [WB]
- 3) Corcelle-Termeau, E., *et al.*, *Autophagy* **12**, 833-849 (2016) [WB]
- 4) Velikkakath, A. K., *et al.*, *Mol. Biol. Cell* **23**, 896-909 (2012) [WB, IC]
- 5) Suzuki, K., *et al.*, *Genes Cell* **12**, 209-218 (2007)
- 6) Suzuki, K., *et al.*, *EMBO J.* **20**, 5971-5981 (2001)
- 7) Strømhaug, P. E., *et al.*, *J. Biol. Chem.* **276**, 42422-42435 (2001)
- 8) Shintani, T., *et al.*, *J. Biol. Chem.* **276**, 30452-30460 (2001)
- 9) Wang, C. W., *et al.*, *J. Biol. Chem.* **276**, 30442-30451 (2001)

For more information, please visit our web site <http://ruo.mbl.co.jp/>

RELATED PRODUCTS

Antibodies

PM036	Anti-LC3 pAb	[WB, IP, IC, IHC, FCM]
M152-3	Anti-LC3 mAb (4E12)	[WB, IP, IC, FCM, EM]
M186-3	Anti-LC3 mAb (8E10)	[WB]
M186-7	Anti-LC3 mAb-HRP-Direct (8E10)	
PD014	Anti-LC3 pAb	[WB]
PM045	Anti-p62 (SQSTM1) pAb	
M162-3	Anti-p62 (SQSTM1) (Human) mAb (5F2)	
M162-A48	Anti-p62 (SQSTM1) (Human) mAb -Alexa Fluor [®] 488 mAb (5F2)	
M162-A59	Anti-p62 (SQSTM1) (Human) mAb -Alexa Fluor [®] 594 (5F2)	
M162-A64	Anti-p62 (SQSTM1) (Human) mAb -Alexa Fluor [®] 647 (5F2)	
PM066	Anti-p62 C-terminal pAb	
D343-3	Anti-Phospho-p62 (SQSTM1) (Ser403) mAb (4F6)	
D344-3	Anti-Phospho-p62 (SQSTM1) (Ser403) mAb (4C8)	
PM074	Anti-Phospho-p62 (SQSTM1) (Ser351) pAb	
M217-3	Anti-Phospho-p62 (SQSTM1) (Ser351) mAb (5D5)	
PD017	Anti-Becclin 1 pAb	
PM037	Anti-GABARAP pAb	
M135-3	Anti-GABARAP mAb (1F4)	
PM038	Anti-GATE-16 pAb	
PD041	Anti-Atg2A pAb	
M133-3	Anti-Atg3 mAb (3E8)	
M134-3	Anti-Atg4B mAb (9H5)	
PM050	Anti-Atg5 pAb	
M153-3	Anti-Atg5 mAb (4D3)	
PM039	Anti-Atg7 (Human) pAb	
PD042	Anti-Atg9A pAb	
M151-3	Anti-Atg10 (Human) mAb (5A7)	
M154-3	Anti-Atg12 (Human) mAb (6E5)	
M183-3	Anti-Atg13 mAb (5G4)	
PD026	Anti-Atg14 pAb	
M184-3	Anti-Atg14 (Human) mAb (4H8)	
PM040	Anti-Atg16L pAb	
M150-3	Anti-Atg16L mAb (1F12)	
M160-3	Anti-UVRAG mAb (1H4)	
PD027	Anti-Rubicon (Human) pAb	
M170-3	Anti-Rubicon (Human) mAb (1H6)	
PM072	Anti-VMP1 pAb	
PM076	Anti-Syntaxin-17 (Human) pAb	
M212-3	Anti-Syntaxin-17 (Human) mAb (2F8)	
PM069	Anti-NRF2 pAb	
M200-3	Anti-NRF2 mAb (1F2)	
M224-3	Anti-KEAP1 mAb (KP1)	
M230-3	Anti-Parkin mAb (Par6)	

WB: Western blotting
IP: Immunoprecipitation
IC: Immunocytochemistry
IHC: Immunohistochemistry
FCM: Flow cytometry
EM: Immuno-electron microscopy

Other related antibodies and kits are also available.
Please visit our website at <http://ruo.mbl.co.jp/>


Kits

8485	Autophagy Ab Sampler Set
8486	Autophagy Watch
CY-7055	CycLex [®] Total p62 ELISA Kit
CY-7056	CycLex [®] Phospho-p62 Ser349 ELISA Kit
CY-7057	CycLex [®] Phospho-p62 Ser403 ELISA Kit
PM036-PN	Positive control for anti-LC3 antibody

SDS-PAGE & Western blotting

- 1) Wash 1×10^7 cells 3 times with PBS and suspend them in 1 mL of Laemmli's sample buffer, then sonicate briefly (up to 20 sec.).
- 2) Centrifuge the tube at $12,000 \times g$ for 5 min. at 4°C and transfer the supernatant to another tube.
- 3) Boil the samples for 3 min. and centrifuge. Load $10 \mu\text{L}$ of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel (7.5% acrylamide) for electrophoresis.
- 4) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at $1 \text{ mA}/\text{cm}^2$ for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 5) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) overnight at 4°C .
- 6) Wash the membrane with PBS-T (0.05% Tween-20 in PBS) (5 min. x 3 times).
- 7) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) as suggested in the **APPLICATIONS** for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 8) Wash the membrane with PBS-T (5 min. x 3 times).
- 9) Incubate the membrane with 1:10,000 Anti-IgG (Rabbit) pAb-HRP (MBL; code no. 458) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 10) Wash the membrane with PBS-T (5 min. x 3 times).
- 11) Wipe excess buffer on the membrane, and then incubate it with appropriate chemiluminescence reagent for 1 min. Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 12) Expose to an X-ray film in a dark room for 1 min. Develop the film as usual. The condition for exposure and development may vary.

(Positive controls for Western blotting; 293T, MEF, PC12 and CHO)




Western blot analysis of Atg2A

Lane 1: 293T
Lane 2: MEF
Lane 3: PC12
Lane 4: CHO
Immunoblotted with PD041

Immunoprecipitation

- 1) Wash 1×10^7 cells 2 times with PBS and resuspend them with 1 mL of ice-cold Lysis buffer [50 mM Tris-HCl (pH 7.5), 150 mM NaCl, 0.05% NP-40] containing appropriate protease inhibitors, then sonicate briefly (up to 20 sec.).
- 2) Centrifuge the tube at 12,000 x g for 5 min. at 4°C and transfer the supernatant to another tube.
- 3) Mix 20 μ L of 50% protein A agarose beads slurry resuspended in 300 μ L of IP buffer [10 mM Tris-HCl (pH 8.0), 500 mM NaCl, 0.1% NP-40] with primary antibody as suggested in the **APPLICATIONS**. Incubate with gentle agitation for 1 hr. at room temperature. (The amount of antibody will depend on the conditions.)
- 4) Centrifuge the tube at 2,500 x g for 10 seconds and discard the supernatant.
- 5) Resuspend the beads with 1 mL of IP buffer.
- 6) Centrifuge the tube at 2,500 x g for 10 seconds and discard the supernatant.
- 7) Repeat steps 4)-6) 2 times.
- 8) Add 300 μ L of cell lysate (prepared sample from step 2)), then incubate with gentle agitation for 1 hr. at room temperature.
- 9) Wash the beads 5 times with 1 mL of Lysis buffer.
- 10) Resuspend the beads in 20 μ L of Laemmli's sample buffer, boil for 3 min. and centrifuge.
- 11) Load 10 μ L of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel (7.5% acrylamide) for electrophoresis.
- 12) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm² for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 13) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) overnight at 4°C.
- 14) Wash the membrane with PBS-T (0.05% Tween-20 in PBS) (5 min. x 3 times).
- 15) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) as suggested in the **APPLICATIONS** for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 16) Wash the membrane with PBS-T (5 min. x 3 times).
- 17) Incubate the membrane with 1:10,000 Anti-IgG (Rabbit) pAb-HRP (MBL; code no. 458) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 18) Wash the membrane with PBS-T (5 min. x 3 times).
- 19) Wipe excess buffer on the membrane, and then incubate it with appropriate chemiluminescence reagent for 1 min.
- 20) Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 21) Expose to an X-ray film in a dark room for 1 min. Develop the film as usual. The condition for exposure and development may vary.

(Positive control for Immunoprecipitation; 293T)



Immunoprecipitation of Atg2A from 293T

Lane 1: IP with Normal Rabbit IgG (PM035)


Lane 2: IP with PD041

Immunoblotted with PD041

Immunocytochemistry

- 1) Spread the cells on a glass slide, then incubate in a CO₂ incubator for one night.
- 2) Remove the culture supernatant by careful aspiration.
- 3) To obtain serum-starved conditions, culture the cells with Hanks' balanced salt solution or DMEM for 2-4 hr. at 37°C.
- 4) Fix the cells by immersing the slide in 4% paraformaldehyde (PFA)/PBS for 10 min. at room temperature (20~25°C).
- 5) Wash the slide 2 times with PBS.
- 6) Immerse the slide in 100 µg/mL of Digitonin in PBS for 10 min. at room temperature.
- 7) Wash the slide 2 times with PBS.
- 8) Add 200 µL of the primary antibody diluted with PBS as suggested in the **APPLICATIONS** onto the cells and incubate for 1 hr. at room temperature. (Optimization of antibody concentration or incubation condition is recommended if necessary.)
- 9) Wash the slide 2 times with PBS.
- 10) Add 200 µL of 1:500 anti-IgG (Rabbit)-Alexa Fluor[®]488 (Invitrogen; code no. A11008) diluted with PBS onto the cells. Incubate for 30 min. at room temperature. Keep out light by aluminum foil.
- 11) Wash the slide 2 times with PBS.
- 12) Promptly add mounting medium onto the slide, then put a cover slip on it.

(Positive control for Immunocytochemistry; starved HeLa)



Immunocytochemical detection of Atg2A in starved HeLa