Ab-Match Universal kit

For technical material or related information, please refer to <https://ruo.mbl.co.jp/bio/product/abmatch/abmatch.html>.

■ **Advantage**

Ab-Match Universal kit comprises of one 96-well microplate, some buffers, some diluents, the substrate solution and the stop solution. Please use Ab-Match Assembly kit in combination with Ab-Match Universal kit in order to build the ELISA system and quantify the target protein in suitable biological samples.

The two kits are sufficient to produce one 96-well microplate.

This product is intended for research use only.

■ **Assay principal**

By using the Ab-Match Universal kit in combination with the Ab-Match Assembly kit for target protein, the sandwich ELISA system is build.

The capture antibody is coated on a 96-well microplate. Standards or samples are added to the microwells, allowing target protein to bind to the capture antibody in proportion to the concentration of target protein in the samples. After wash unbound target protein, detection antibody is added and attaches to the target protein bound to coated antibody. After washing away all unbound detection antibody, initiate the substrate reaction. After substrate reaction, add the stop solution to the microwell plate. Then the relative amount of target protein bound in each well can be determined from the yellow color using a microplate reader. The concentration of target protein in each sample can be obtained by comparing the optical density (OD) of the sample to the OD of the standards on the plate.

■ **Materials provided** 96 wells

Ab-Match Universal kit contains:

Name	Materials	Quantity (for 1 kit)
Microwell strips	Non coated microwell strips	8-well × 12 strips
Coating Buffer	Carbonate buffer solution (ready-to-use)	12 mL × 1 bottle
Blocking Agent	Contains BSA and sucrose (ready-to-use)	24 mL × 1 bottle
Sample Diluent	Buffer mostly for diluting samples (ready-to-use) Contains BSA, Tween 20 and HAMA-Blocker	50 mL × 1 bottle
Wash Concentrate	Buffer for washing microwells (20×) Contains Tween 20	50 mL × 1 bottle
SA-HRP Diluent	Buffer mostly for diluting SA-HRP (ready-to-use) Contains BSA	15 mL × 1 bottle
Substrate Solution	TMB/H ₂ O ₂ solution (ready-to-use)	20 mL × 1 bottle
Stop Solution	0.25 mol/L sulfuric acid (ready-to-use)	20 mL × 1 bottle

Note: Kit components must be stored at 2-8°C. Microwell strips can be stored at room temperature.

Kit components must be brought to room temperature (20-30°C) prior to use.

■ **Precautions**

- Do not use reagents beyond the stated expiration dates.
- A standard curve must be run with each assay.
- Operation for dispensing and diluting should be precisely done.
- Avoid contact of **Substrate Solution** and **Stop Solution** with skin or eyes. If contacted, wash away with plenty of water.
- **Substrate Solution** is easily oxidized with metal ions. Use disposable new instruments and disposable pipettes for all handling of the substrate solution NEVER return **Substrate Solution** to the bottle.
- Serum samples may be infectious. Instruments used in this test should be disposed after use or treated as follows:
 - Soak in 2% glutaraldehyde solution (final concentration) for more than one hour.
 - Soak in 0.5% sodium hypochlorite solution (available chloric: approximately 5,000 ppm.) for more than one hour.
 - Autoclave at 121°C for more than 20 minutes.
- The incubation times indicated do not allow the incubation to complete. Frequent moving of the plate during standing, or vibration from instruments, may cause a shaking effect to the reaction solution, causing the reaction to progress faster than usual and giving higher color development.


■ **Procedure**

Prepare and use the Ab-Match Assembly kit reagent's standard and antibody solution with Ab-Match Universal kit's buffer solution. Ab-Match Universal kit contains the Microwell strips, Coating Buffer, Sample Diluent, Wash Concentrate, SA-HRP Diluent, Substrate Solution, and Stop Solution. These are required to prepare one 96-well ELISA using one Ab-Match Assembly kit.

Kit components must be brought to room temperature (20-30°C) prior to use.

The details of assay procedure vary from Ab-Match Assembly kit to kit. Please read data sheet of each Ab-Match Assembly kit thoroughly.

◆ **Outline**

* Please read data sheet of each Ab-Match Assembly kit thoroughly.

◆ **Preparation of Reagents**

For reagents highlighted in a box (), use **Ab-Match Universal kit** components. Each vial should be centrifuged prior to use in order to collect the reagent that is on the wall and cap of each vial.

1. Wash Concentrate

Prepare **Wash solution** by diluting **Wash Concentrate**, 1:20 with the deionized water prior to use (ex. Add 50 mL of **Wash Concentrate** to 950 mL of deionized water). After preparation, it is stable for 2 weeks at 2-8 °C.

2. standard (component of Ab-Match Assembly kit)

For preparing **Conditioned standard**, **standard** should be diluted 1:10 with **Sample Diluent** or culture medium. With the 1:10 dilution as top **Conditioned standard**, 5 to 7 points should be made with **Sample Diluent** or culture medium by 2-fold serial dilution. Use the same diluent without any standard added as blank (zero conc.).

“**standard**” must be stored below -20°C. To avoid repeated freezing and thawing, prepare aliquots.

3. capture antibody (component of Ab-Match Assembly kit)

The capture antibody should be diluted with Coating Buffer.

Dilution must be used immediately and cannot be stored for future use.

Please read data sheet of each Ab-Match Assembly kit thoroughly.

4. detection antibody (component of Ab-Match Assembly kit)

The Biotin conjugated detection antibody should be diluted with Sample Diluent.

The HRP conjugated detection antibody should be diluted with SA-HRP Diluent.

Dilution must be used immediately and cannot be stored for future use.

Please read data sheet of each Ab-Match Assembly kit thoroughly.

5. Streptavidin conj. peroxidase (component of some Ab-Match Assembly kits)

In case **SA-HRP solution** is needed, prepare **SA-HRP solution** by diluting Streptavidin conj. peroxidase (SA-HRP) with SA-HRP Diluent.

Dilution must be used immediately and cannot be stored for future use.

Please read data sheet of each Ab-Match Assembly kit thoroughly.

◆ Microplate Coating with capture antibody

For reagents highlighted in a box (), use components of the **Ab-Match Universal kit**.

STEP 1. (Coating)

1) Dilute capture antibody with Coating Buffer.

2) Soon after mixing by repeated inversion, dispense 100 µL to each microwell with multichannel pipette and let stand overnight at 2-8°C with seal (or other cover to avoid evaporation).

* Use a new conical 15mL tube to mix the antibody solution.

* Please confirm the optimal dilution ratio with the data sheet of the Ab-Match Assembly kit.

STEP 2. (Washing)

Wash the wells with saline (saline: NaCl 9.0g / 1,000 mL) 2 times.

Washing Method

Discard the liquid out of the wells. Be careful not to detach strips from the microcup frame.

Gently pour wash fluid into the empty wells.

Repeat the discarding and refilling steps.

Finally, tap the plate several times on paper towel to remove any contents.

* If an autowasher is used, optimal washing times vary depending on the instrument used and its setting.

* Wash as quickly as possible, NOT let wells dry up.

STEP 3. (Blocking)

1) Add 200 µL of Blocking Agent to each well and incubate for one hour at room temperature.

2) Dump out the contents of the wells over sink before use.

* Discarding of antibody and addition of blocking should be done as soon as possible as to not to let the cup become dry.

* After completion of the blocking step, the coated plate can be stored for a long period of time when properly dried. After discarding blocking agent, apply fan or other drying method and leave at room temperature for 3 hours to overnight to ensure proper drying. Store the plate at 2-8°C, under strictly controlled moisture conditions. The color development of dried plates may be decreased when compared with plates that were never dried.

◆ **Assay procedure**

Duplicate assay is recommended. A standard curve must be run with each assay.

STEP 1. (Sample incubation)

1) Dilute standard 1:10 with Sample Diluent or culture medium.

For preparing dilution series, with the 1:10 dilution as top **Conditioned standard**, 5 to 7 points should be made with Sample Diluent or culture medium by 2-fold serial dilution. Then use buffer solution used for preparation as Blank.

2) Dilute samples with optimal dilution in Sample Diluent or culture medium.

* When assayed with culture supernatant as sample, culture medium should be used to dilute the standard.

* Sample Diluent contains HAMA-Blocker to block the effect of human anti- mouse antibodies (HAMA) that may be present in human serum.

3) Add 100 µL of each **Conditioned standard** points and samples to coated well.

Incubate for 1 hour at room temperature (primary reaction)

* The Antigen-antibody reaction starts on addition. Addition should be completed as quickly as possible. It is recommended that standard and sample are diluted on a separate microplate in advance, then added to the antibody coated plate with a multichannel pipette.

* When adding solution to the microplate, avoid touching the inner wall of microcup with the pipet tip. This technique avoids non-specific reaction.

STEP 2. (Washing)

Wash the wells 4 times with **Wash solution** following the Washing Method.

For Washing Method, refer to “◆ **Microplate Coating with capture antibody**, STEP 2”.

STEP 3. (Detection antibody incubation)

1) Add 100 µL of diluted detection antibody to each well.

2) Incubate for 1 hour at room temperature.

* Avoid touching the inner wall of the wells with the pipette tip. Otherwise, it could cause of non-specific reaction.

* Please confirm the optimal dilution ratio with the data sheet of the Ab-Match Assembly kit.

STEP 4. (Washing)

Wash the wells 4 times with **Wash solution** following the Washing Method.

STEP 5. (SA-HRP incubation)

In case that the detection antibody is conjugated with HRP, STEP 5 and 6 are not needed.

1) Add 100 µL of diluted SA-HRP solution to each well.

2) Incubate for 30 minutes at room temperature.

* Avoid touching the inner wall of the wells with the pipette tip. Otherwise, it could cause of non-specific reaction.

STEP 6. (Washing)

Wash the wells 4 times with **Wash solution** following the Washing Method.

STEP 7. (Substrate incubation)

Add 100 µL of Substrate Solution to each well. Incubate for 30 minutes at room temperature.

STEP 8. (Stopping reaction)

Add 100 µL of Stop Solution to each well.

CODE No. 5310

◆ **Reading**

Using microplate reader, read the absorbance of each well at 450 nm and the reference at 620 nm.

* Read absorbance within 30 minutes after reaction stop.

* Ensure that the back of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading.

◆ **Calculation of result**

Calculate the mean absorbance value of each standard and plot against log standard concentration and connect the points the best fitting straight line. The concentration of the samples then can be read from this standard curve. Alternatively a suitable computer and curve-fitting program can be used. The concentration read from the standard curve must be multiplied by the dilution factor.

* If a sample's OD is out of range for the standard curve, the assay should be repeated with a higher sample dilution.

ODs should always remain below 2.0 in order to remain in the dynamic range of the detection system.

* When computer software is used, logistic, 3 (4) -para-logit-log, or Spline may be used as calculation.

* A concentration conversion spreadsheet using 4-para-Logistic regression with Excel 2000 (Microsoft) is available

Calculation Program Operating Manual: https://ruo.mbl.co.jp/bio/product/abmatch/images/abmatch/manual_e.pdf

Calculation Program :https://ruo.mbl.co.jp/bio/product/abmatch/images/abmatch/logistic_e_v2.xls

■ **Storage**

All kit components must be stored at 2~8°C.

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Ab-Match Universal kit:

Name	Materials	Quantity (for 1 kit)
Microwell strips	Non coated microwell strips	8-well × 12 strips
Coating Buffer	Carbonate buffer solution (ready-to-use)	12 mL × 1 bottle
Blocking Agent	Contains BSA and sucrose (ready-to-use)	24 mL × 1 bottle
Sample Diluent	Buffer for diluting samples (ready-to-use) Contains BSA, Tween 20 and HAMA-Blocker	50 mL × 1 bottle
Wash Concentrate	Buffer for washing microwells (20×) Contains Tween 20	50 mL × 1 bottle
SA-HRP Diluent	Buffer for diluting SA-HRP (ready-to-use) Contains BSA	15 mL × 1 bottle
Substrate Solution	TMB/H ₂ O ₂ solution (ready-to-use)	20 mL × 1 bottle
Stop Solution	0.25 mol/L sulfuric acid (ready-to-use)	20 mL × 1 bottle

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MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.