

**For Research Use Only.
Not for use in diagnostic procedures.**

Afamin/Wnt3a CM (Conditioned Medium)

CODE No.	J2-001
FORM	Liquid
QUANTITY	10 mL
DESCRIPTION	<p>Conditioned medium containing Wnt3a protein and Afamin protein.</p> <p>The mouse Wnt3a and the human Afamin were co-expressed in Expi293F™ Cells*. For the co-expression of Wnt3a and Afamin, expression plasmids for N-terminal PA tagged mouse wnt3a and N-terminal TARGET tagged human afamin were co-transfected using the Expi293 Expression System*.</p> <p>Conditioned medium was prepared using serum-free media (Advanced DMEM/F-12* containing penicillin, streptomycin, HEPES and GlutaMAX™ Supplement*). Then, it was sterilized by filtration of 0.2 µm-filter.</p> <p>FBS was used for pre-culture of Afamin/Wnt3a expression Expi293F™ Cells.</p> <p>*: Thermo Fisher Scientific's products.</p>
ENDOTOXIN LEVEL	< 1.0 EU/mL by the LAL assay.
STORAGE	This conditioned medium is stable for 18 months from the date of manufacture when stored at -20°C or below.
DIRECTION for USE	<p>As a canonical Wnt signal agonist.</p> <p>For intestinal organoids culture, add to the culture media at a final concentration of 10%. Since the optimal condition may vary depending on samples, each user should determine the optimal volume of addition in advance.</p>

BACKGROUND

Wnt signaling is known to be involved in early development, maintenance and regeneration of stem cells, and in cancer formation. Wnt signaling has also been found to play an important role in the growth and maintenance of these processes. In particular, Wnt3a has been revealed to be an essential niche component for maintaining the proliferation of Lgr5-positive stem cells in intestinal epithelial cells and is used for the production of various digestive organoids such as the small intestine, large intestine, stomach, pancreas and liver. Although Wnt3a has been conventionally used for the culture of gut organoids, it is a fat-soluble protein, so it forms aggregates in serum-free medium and can not exert its activity sufficiently. In 2016, Mihara *et al.* found that high Wnt3a activity can be maintained by forming a complex with Wnt3a by Afamin, which is one of the components of serum. In addition, by using Afamin and Wnt3a complex for organoid culture, long-term culture of organoid becomes possible. This new medium will result in optimal success for your organoid experiments.

REFERENCES

- 1) E. Mihara, *et al.*, Active and water-soluble form of lipidated Wnt protein is maintained by a serum glycoprotein afamin/ α -albumin., *eLife* **5** (2016) [PMID: [26902720](#)]
- 2) K. Nanki, *et al.*, Divergent routes toward Wnt and R-spondin niche independency during human gastric carcinogenesis., *Cell* **174** (2018) [PMID: [30096312](#)]
- 3) S. Sugimoto, *et al.*, Reconstruction of the human colon epithelium in vivo., *Cell Stem Cell* **22** (2018) [PMID: [29290616](#)]
- 4) T. Seino, *et al.*, Human pancreatic tumor organoids reveal loss of stem cell niche factor dependence during disease progression., *Cell Stem Cell* **22** (2018) [PMID: [29337182](#)]
- 5) S. Sugimoto, *et al.*, Organoid Derivation and Orthotopic Xenotransplantation for Studying Human Intestinal Stem Cell Dynamics., *Methods Mol Biol.* **2171** (2020) [PMID: [32705652](#)]
- 6) K. Nanki, *et al.*, Somatic inflammatory gene mutations in human ulcerative colitis epithelium., *Nature* **577** (2020) [PMID: [31853059](#)]
- 7) N. Sasaki, *et al.*, Development of a Scalable Coculture System for Gut Anaerobes and Human Colon Epithelium. , *Gastroenterology* **159** (2020) [PMID: [32199883](#)]
- 8) S. Mae, *et al.*, Expansion of Human iPSC-Derived Ureteric Bud Organoids with Repeated Branching Potential. ,*Cell Reports* **32** (2020)[PMID: [32726627](#)]
- 9) Y. Nanki, *et al.*, Patient-derived ovarian cancer organoids capture the genomic profiles of primary tumours applicable for drug sensitivity and resistance testing. ,*Scientific Reports* **28** (2020) [PMID: [32724113](#)]
- 10) K. Miyabayashi, *et al.*, Intraductal transplantation models of human pancreatic ductal adenocarcinoma reveal progressive transition of molecular subtypes. , *Cancer Discov.* **10** (2020)[PMID: [32703770](#)]

For more information, please visit our web site <https://ruo.mbl.co.jp/>.

When culturing organoids, or stem cells, or other tissues, if you are to use this product in combination with other factor or factors (hereunder factors), a third party may have a patent on the use or other application of the factors concerned. Regarding to this product, we do not offer any non-infringement warranty when used or otherwise applied in combination with other factors. Therefore, if you intend to use this product in combination with other factors, please check with your organization's division responsible for intellectual property rights or your research agency before using this product.