For Research Use Only. Not for use in diagnostic procedures.

# **DDDDK-tagged Protein**

## **Magnetic PURIFICATION KIT**

(mAb. clone FLA-1)

CODE No. 3343

*PURIFICATION to Maintain Protein Activity  
from eukaryote cell lysate and culture supernatant*

## Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more.

The method of purification is one of the important keys for maintaining protein structure and function. The DDDDK-tagged Protein Magnetic PURIFICATION KIT is designed for the isolation of DDDDK-tagged protein from cell lysates and culture supernatant. Severe conditions such as acidic or alkaline elution can denature protein structure while maintaining a neutral pH can preserve protein activity and conformation. MBL has developed the Anti-DDDDK-tag Magnetic Beads to quickly and efficiently purify DDDDK-tagged proteins at neutral pH to maintain protein activity and native conformation. The DDDDK-tag peptide competitively elutes the DDDDK tagged protein from Anti-DDDDK-tag Magnetic Beads while minimally interfering with protein function.

The affinity of the anti-DDDDK-tag antibody is also very important. Very high affinity antibody cannot be dissociated from the DDDDK-tagged protein, and low affinity antibody is not sufficient for binding. MBL has optimized the affinity of the anti-DDDDK-tag antibody to permit high yields and efficient purification of functional DDDDK-tagged proteins for biochemical studies. This kit is simple and fast, because all procedures of this kit have been optimized.

## Kit Components.

### 1. Anti-DDDDK-tag Magnetic Beads 1 mL x 2 vials

15 mg magnetic beads in 1 mL PBS/0.1% BSA with 0.09% sodium azide as preservative

\*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.


### 2. Elution Peptide 1 mg x 2 vials (lyophilized)

DDDDK-tag peptide (DYKDDDDK), 1 mg in 1 mL PBS after reconstitution

### 3. Wash Concentrate 5 mL x 2 vials

20x concentrate

## Procedure Summary

## Product Capacity

The purification capacity of the Anti-DDDDK-tag Magnetic Beads varies depending upon the DDDDK tagged protein. For examples, 100  $\mu$ L of Anti-DDDDK-tag Magnetic Beads (100  $\mu$ L slurry) can be purified 2.0  $\mu$ g of a DDDDK-tagged protein (32 kDa).

## **Materials Required but not Provided**

1. Microcentrifuge capable of 15,000 x g
2. Sampling tube (1.5 mL)
3. End-over-end rotator
4. Magnetic Rack (MBL; code 3190)
5. Lysis buffer

Suitable Lysis buffer varies with cell type.

**Note:** see Additional Information

### Homemade Lysis buffer

|           |   |
|-----------|---|
| 20-50 mM  | Tris-HCl (pH 7.5) or HEPES-KOH (pH 7.5)   |
| 50-250 mM | NaCl  |
| 5 mM      | EDTA  |
| 1%        | NP-40 or Triton X-100   |
|           | if necessary add Protease Inhibitor Cocktail<br>(e.g. SIGMA; code P8340, PIERCE; code 78415). |

## **Storage**

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

## **Protocols**

The following protocol is for the isolation of DDDDK-tagged proteins produced in a 100-mm cell culture dish. The expression level of the DDDDK-tagged protein may vary. If necessary, adjust the volume of Anti-DDDDK-tag Magnetic Beads and Elution Peptide Solution proportionally.

### **Material Preparation**

#### 1. Wash Solution

Dilute Wash Concentrate with 20 times its volume of distilled water. (e.g. Dilute 0.5 mL of Wash Concentrate with 9.5 mL of distilled water.) For each purification, prepare 5 mL of Wash Solution.

#### 2. Elution Peptide Solution

Reconstitute the Elution Peptide with 1 mL of distilled water. If you want to store the reconstituted Elution Peptide, prepare appropriate aliquots (e.g. 45 µL x 20 tubes) and store at -20°C. Repeated freezing and thawing is not recommended.

## A. Purification from mammalian cultured cell lysate

(DDDDK-tagged protein is not secreted from the cells)

### (Lysis of Mammalian Cells)

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 0.5 mL of Lysis buffer to the cell pellet and vortex.
6. Incubate the sample for 15 minutes on ice.
7. Remove cell debris by centrifugation at 15,000 x g for 5 minutes at 4°C.

### (Purification of DDDDK-tagged Protein)

8. Transfer the 0.5 mL of cell lysate to the 1.5 mL sampling tube.
9. Resuspend Anti-DDDDK-tag Magnetic Beads by vortexing for a few seconds immediately before dispensing.
10. Dispense 50 µL of the Anti-DDDDK-tag Magnetic Beads suspension into the sampling tube prepared in step 8. If necessary, please use 100 µL of the Anti-DDDDK-tag Magnetic Beads suspension to increase protein yield.
11. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
12. Place the tube on Magnetic Rack (MBL; code no.3190) for 10 seconds. And remove the supernatant.
13. Add 1 mL of Wash Solution to the tube. Resuspend the magnetic bead pellet by vortexing for a few seconds.
14. Place the tube on the Magnetic Rack for 10 seconds. And remove the supernatant.
15. Repeat steps 13-14 three additional times.
16. Add 20 µL Elution Peptide Solution to the tube. Vortex the tube for a few seconds. Incubate for 5 minutes at room temperature.
17. Place the tube on the Magnetic Rack for 10 seconds. Transfer the supernatant to another tube.
18. Repeat steps 16-17 again. The two eluates may be pooled in one sampling tube.

## B. Purification from culture supernatant

(DDDDK-tagged protein is secreted into the culture supernatant)

1. Collect the culture supernatant from the cell culture dish into 15 mL centrifuge tube.
2. Centrifuge at 400 x g for 5 minutes to remove cell debris.
3. Transfer the supernatant to new 15 mL centrifuge tube.
4. Resuspend Anti-DDDDK-tag Magnetic Beads by vortexing for a few seconds immediately before dispensing.
5. Dispense 50 µL of the Anti-DDDDK-tag Magnetic Beads suspension into the 15 mL centrifuge tube with culture supernatant. If necessary, please use 100 µL of the Anti-DDDDK-tag Magnetic Beads suspension to increase protein yield.
6. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
7. Centrifuge at 400 x g for 5 minutes using swinging bucket rotor. Carefully remove the supernatant, but leave 100-400 µL supernatant above the magnetic bead pellet.
8. Resuspend the magnetic bead pellet in the 100-400 µL supernatant by vortexing for a few seconds.

- Transfer the suspension to 1.5 mL sampling tube.
9. Place the tube on Magnetic Rack (MBL; code no. 3190) for 10 seconds. Remove the supernatant.
  10. Add 1 mL of Wash Solution to the tube. Resuspend the magnetic bead pellet by vortexing for a few seconds.
  11. Place the tube on the Magnetic Rack for 10 seconds. Remove the supernatant.
  12. Repeat steps 10-11 three additional times.
  13. Add 20  $\mu$ L Elution Peptide Solution to the tube. Vortex the tube for a few seconds. Incubate for 5 minutes at room temperature.
  14. Place the tube on the Magnetic Rack for 10 seconds. Transfer the supernatant to another tube.
  15. Repeat steps 13-14 again. The two eluates may be pooled in one sampling tube.


(Note: Steps 14-18 of protocol A are identical to steps 11-15 of protocol B.)

### **Related Products:**

|          |   |                                    |
|----------|---|------------------------------------|
| 3325     | DDDDK-tagged Protein Purification Kit                           | 20 purifications                   |
| 3326     | DDDDK-tagged Protein Purification Gel with Elution Peptide      |                                    |
|          |   | Gel: 1 mL x 1, Peptide: 1 mg x 5   |
| 3328     | DDDDK-tagged Protein Purification Gel                           | 5 mL x 1                           |
| 3325-205 | DDDDK-tag peptide   | 1 mg x 5                           |
| 3305     | c-Myc-tagged Protein Mild Purification Kit ver.2                | 20 purifications                   |
| 3305A    | c-Myc-tagged Protein Mild Purification Kit ver.2                | 2 purifications                    |
| 3306     | c-Myc-tagged Protein Mild Purification Gel with Elution Peptide |                                    |
|          |   | Gel: 1 mL x 1, Peptide: 1 mg x 5   |
| 3307     | c-Myc-tagged Protein Mild Purification Gel with Elution Peptide |                                    |
|          |   | Gel: 1 mL x 5, Peptide: 1 mg x 5   |
| 3300-205 | c-Myc-tag peptide   | 1 mg x 5                           |
| 3340     | c-Myc-tagged Protein Magnetic Purification Kit                  |                                    |
|          |   | Beads: 1 mL x 2, Peptide: 1 mg x 2 |
| 3310     | His-tagged Protein Purification Kit                             | 20 purifications                   |
| 3311     | His-tagged Protein Purification Gel with Elution Peptide        |                                    |
|          |   | Gel: 1 mL x 1, Peptide: 2 mg x 5   |
| 3310-205 | His-tag peptide   | 2 mg x 5                           |
| 3317     | V5-tagged Protein Purification Kit Ver.2                        | 20 purifications                   |
| 3318     | V5-tagged Protein Purification Gel Ver.2                        | 1 mL x 1                           |
| 3315-205 | V5-tag peptide  |                                    |
| 3341     | V5-tagged Protein Magnetic Purification Kit                     |                                    |
|          |   | Beads: 1 mL x 2, Peptide: 2 mg x 2 |
| 3320     | HA-tagged Protein Purification Kit                              | 20 purifications                   |
| 3321     | HA-tagged Protein Purification Gel                              | 1 mL                               |
| 3320-205 | HA-tag peptide  | 2 mg x 5                           |
| 3342     | HA-tagged Protein Magnetic Purification Kit                     |                                    |
|          |   | Beads: 1 mL x 2, Peptide: 2 mg x 2 |

## Example of Purification Result

### Purification of DDDDK-tagged $\beta$ -galactosidase

Lane 1: Input

Lane 2: Eluate (purified using DDDDK-tagged Protein Magnetic PURIFICATION KIT)

Human embryonic kidney (HEK293T) cells were transfected with pcDNA-DDDDK-tagged  $\beta$ -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (0.5 mL/100-mm dish) and purified according to the preceding protocol A.

### **Additional Information**

Several reagents were examined whether or not they were suitable for use with the DDDDK-tagged Protein Magnetic PURIFICATION KIT. The results are listed below.

#### **Chaotropic agents**

|                      |            |            |
|----------------------|------------|------------|
| <b>Urea</b>          | <b>1 M</b> | <b>Yes</b> |
| <b>Guanidine-HCl</b> | <b>1 M</b> | <b>No</b>  |

#### **Reducing agents**

|                          |              |            |
|--------------------------|--------------|------------|
| <b>DTT</b>               | <b>10 mM</b> | <b>Yes</b> |
| <b>2-Mercaptoethanol</b> | <b>10 mM</b> | <b>Yes</b> |

#### **Surfactants**

|                     |   |              |            |
|---------------------|---|--------------|------------|
| <b>Nonionic</b>     | <b>Tween-20</b>                               | <b>5%</b>    | <b>Yes</b> |
|                     | <b>TritonX-100</b>                            | <b>5%</b>    | <b>Yes</b> |
|                     | <b>NP-40</b>                                  | <b>5%</b>    | <b>Yes</b> |
|                     | <b>Digitonin</b>                              | <b>1%</b>    | <b>Yes</b> |
|                     | <b>n-Octyl-<math>\beta</math>-D-gulcoside</b> | <b>1%</b>    | <b>Yes</b> |
| <b>Zwitterionic</b> | <b>CHAPS</b>                                  | <b>1%</b>    | <b>Yes</b> |
|                     | <b>CHAPSO</b>                                 | <b>1%</b>    | <b>Yes</b> |
| <b>Anionic</b>      | <b>SDS</b>                                    | <b>0.05%</b> | <b>Yes</b> |
|                     | <b>Sodium Deoxycholate</b>                    | <b>0.5%</b>  | <b>Yes</b> |

#### **Others**

|                 |              |            |
|-----------------|--------------|------------|
| <b>NaCl</b>     | <b>1 M</b>   | <b>Yes</b> |
| <b>Glycerol</b> | <b>10%</b>   | <b>Yes</b> |
| <b>EDTA</b>     | <b>10 mM</b> | <b>Yes</b> |

The “Yes” indicates the reagents can be used in the Lysis buffer for this kit up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this kit at the indicated concentration.

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
|          |   |                                    |
|----------|---|------------------------------------|
| 3325     | DDDDK-tagged Protein Purification Kit                           | 20 purifications                   |
| 3326     | DDDDK-tagged Protein Purification Gel with Elution Peptide      | Gel: 1 mL x 1, Peptide: 1 mg x 5   |
| 3328     | DDDDK-tagged Protein Purification Gel                           | 5 mL x 1                           |
| 3325-205 | DDDDK-tag peptide   | 1 mg x 5                           |
| 3305     | c-Myc-tagged Protein Mild Purification Kit ver.2                | 20 purifications                   |
| 3305A    | c-Myc-tagged Protein Mild Purification Kit ver.2                | 2 purifications                    |
| 3306     | c-Myc-tagged Protein Mild Purification Gel with Elution Peptide | Gel: 1 mL x 1, Peptide: 1 mg x 5   |
| 3307     | c-Myc-tagged Protein Mild Purification Gel with Elution Peptide | Gel: 1 mL x 5, Peptide: 1 mg x 5   |
| 3300-205 | c-Myc-tag peptide   | 1 mg x 5                           |
| 3340     | c-Myc-tagged Protein Magnetic Purification Kit                  | Beads: 1 mL x 2, Peptide: 1 mg x 2 |
| 3310     | His-tagged Protein Purification Kit                             | 20 purifications                   |
| 3311     | His-tagged Protein Purification Gel with Elution Peptide        | Gel: 1 mL x 1, Peptide: 2 mg x 5   |
| 3310-205 | His-tag peptide   | 2 mg x 5                           |
| 3317     | V5-tagged Protein Purification Kit Ver.2                        | 20 purifications                   |
| 3318     | V5-tagged Protein Purification Gel Ver.2                        | 1 mL x 1                           |
| 3315-205 | V5-tag peptide  |                                    |
| 3341     | V5-tagged Protein Magnetic Purification Kit                     | Beads: 1 mL x 2, Peptide: 2 mg x 2 |
| 3320     | HA-tagged Protein Purification Kit                              | 20 purifications                   |
| 3321     | HA-tagged Protein Purification Gel                              | 1 mL                               |
| 3320-205 | HA-tag peptide  | 2 mg x 5                           |
| 3342     | HA-tagged Protein Magnetic Purification Kit                     | Beads: 1 mL x 2, Peptide: 2 mg x 2 |

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Lane1: Input

Lane2: Eluate (purified using HA-tagged Protein Magnetic PURIFICATION KIT)

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