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**For Research Use Only.  
Not for use in diagnostic procedures.**

**MBL**

## Anti-RFP HRP-Direct

Code No.	Quantity	Antibody type
PM005-7	100 µL	Rabbit Polyclonal

### HRP-Direct Series

HRP-Direct is a series of HRP conjugated primary antibodies developed by MBL.

HRP-Direct Series products don't need secondary antibodies. That brings the following advantages:

① Total incubation time is cut in half.

The reaction with the secondary antibody becomes unnecessary. In the shaking method, reaction time is only 30 minutes. Spend your evening doing more important things.

② Clear result

No more heavy and light chain bands in Immunoprecipitation. The HRP conjugated primary antibody will not detect your precipitating antibody. Clear result helps your research.

Usually, the secondary antibody amplifies the signal of an unconjugated primary antibody. Therefore, it is thought to be difficult to obtain a strong signal with a directly conjugated antibody. To overcome this hurdle, MBL has improved the HRP conjugation method. MBL's HRP-Direct Series yield strong signal with minimum background. Give it a try - the HRP-Direct Series will not disappoint you.


### DETECTED ANTIGEN

RFP (Red Fluorescent Protein) , RFP Variants.

\*Our RFP polyclonal antibody has been confirmed to react with the following RFP Variants.

DsRed, mRFP1, mCherry, mOrange, mPlum, mStrawberry, TdTomato

This antibody does not react with GFP (Green Fluorescent Protein), which distinguish RFP from GFP.

## REAGENT

PBS/Preservative/Stabilizer

## STORAGE

This antibody is stable for one year from the date of purchase when stored at 4°C.

## INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.


## APPLICATION

Western blotting

Method	Reaction times	Dilution	Detaction
Sitting	1 hour	1,000~4,000	chemiluminescence
Shaking	30minutes	1,000~4,000	chemiluminescence

\* Detailed procedure is provided in the following protocols.

Comparison of Sitting method and Shaking method

Immunocytochemistry; Not tested

Immunohistochemistry; Not tested

## Immunoprecipitation (for clear results)

There are two causes of non-specific bands at Western blotting that follows Immunoprecipitation.

### Cause1: Reaction of the secondary antibody with the antibody used at IP. (Fig.3; lane 1)

#### Solution

MBL HRP-Direct series (HRP conjugated primary antibodies products series) does not detect the precipitating antibody (Fig. 3; lanes 2\*, 3, and 4). \*The nonspecific band in lane 2 represents protein A/G. Please refer to **Cause 2**.


### Cause 2: Reaction of protein A/G with the primary and secondary antibodies used at WB. (Fig.3; lane 2)

Boiling the Immunoprecipitated sample with Laemmli's sample buffer; this causes dissociation of protein A/G from its carrier. Therefore, the primary and secondary antibodies bind to protein A/G at western blotting.

#### Solution

- ① Add Clear Back for IP Western (MBL; code no. MTG-002) in 1:100 ratio to the diluted HRP-Direct. (Fig.3; lane 3)
- ② Use an Agarose-conjugated antibody for Immunoprecipitation. (Fig.3; lane 4)

Figure3

## PROTOCOL Western Blotting

### 1. SDS-PAGE and Western Blotting

Perform SDS-PAGE electrophoresis on the protein samples and transfer the protein to a PVDF membrane according to standard techniques.

### 2. Blocking

Block the membrane with 5% skimmed milk in PBS (or TBS). Incubate for 1 hour with shaking at room temperature or over night at 4°C.


3. Incubate the membrane with HRP-Direct antibody as follows:

① Sitting Method (Saving type: amount of antibody)

Dilute the HRP-Direct antibody with 1% skimmed milk in PBS (or TBS). Use at least 30  $\mu\text{L}/\text{cm}^2$  membrane. Antibody dilution ratio is suggested in the APPLICATION. Drain the excess blocking buffer from the membranes and place them, protein side up, on a wrap film or other suitable clean surface. Pipette the HRP-Direct solution on to the membrane. Cover the membrane with plastic case to avoid drying. Incubate for 1 hour at room temperature.

② Shaking Method (Saving type: reaction time)

Dilute the HRP-Direct antibody with 1% skimmed milk in PBS (or TBS). Use at least 10mL. Antibody dilution ratio is suggested in the APPLICATION. Transfer membrane to a tray containing HRP-Direct antibody. Incubate for 30 minutes with shaking at room temperature.

4. Washing

Wash the membrane in PBS (or TBS) with 0.05% Tween-20 three times for 5 minutes each.

5. Detection

Detect with a chemiluminescence substrate according to the manufacturer's instructions.

100514-1.1